

AD A 120835

12

REPORT DRXTH-TE-CR-82179

SURFACE SAMPLING TECHNIQUES

Bruce E. Goodwin
James R. Aronson
Robert P. O'Neil
Margaret A. Randel
Emmett M. Smith

**ARTHUR D. LITTLE, INC.
CAMBRIDGE, MA 02140**

SEPTEMBER, 1982

**FINAL REPORT
Volume II
Certification Testing Data**

**Distribution Unlimited
Cleared for Public Release**

prepared for

U.S. Army Toxic and Hazardous Materials Agency,
Aberdeen Proving Ground, Maryland 21010

**DTIC
ELECTE
OCT 29 1982
S H**

Best Available Copy

copy

DISCLAIMER

The views, opinions and/or findings contained in this report are those of the authors and should not be construed as an official Department of the Army position, policy or decision, unless so designated by other documents.

REPORT DRXTH-TE-CR-82179

SURFACE SAMPLING TECHNIQUES

Bruce E. Goodwin
James R. Aronson
Robert P. O'Neil
Margaret A. Randel
Emmett M. Smith

***ARTHUR D. LITTLE, INC.
CAMBRIDGE, MA 02140***

SEPTEMBER, 1982

**FINAL REPORT
Volume II
Certification Testing Data**

**Distribution Unlimited
Cleared for Public Release**

prepared for

**U.S. Army Toxic and Hazardous Materials Agency,
Aberdeen Proving Ground, Maryland 21010**



TABLE OF CONTENTS :

	<u>Page</u>
Abbreviations	ii
List of Tables	iii
List of Figures	xi
> Semiquantitative Certification Data;	5
Quantitative Method for the Determination of DNP, RDX, TNB, 2,4-DNT, TNT, Tetryl and DPA on Surfaces;	
Quantitative Method for the Determination of 2,6-DNT and NG on Surfaces;	
Quantitative Method for the Determination of PETN on Surfaces; <i>o.v. 6</i>	
Quantitative Concentration Data.	



Accession For	
NTIS GRA&I	<input checked="" type="checkbox"/>
DTIC TAB	<input type="checkbox"/>
Unannounced	<input type="checkbox"/>
Justification	
By	
Distribution/	
Availability Codes	
Avail and/or	
Dist Special	
A	

ABBREVIATIONS

<u>Abbreviation</u>	<u>Chemical Name</u>
DNP	2,4-Dinitrophenol
RDX	Cyclotrimethylenetrinitramine
TNB	1,3,5-Trinitrobenzne
2,4-DNT	2,4-Dinitrotoluene
2,4,6-TNT	2,4,6-Trinitrotoluene
Tetryl	2,4,6-Trinitrophenylmethylnitramine
DPA	Diphenylamine
2,6-DNT	2,6-Dinitrotoluene
NG	Nitroglycerine
PETN	Pentaerythritetetranitrate

LIST OF TABLES

<u>TABLE NO.</u>		<u>PAGE</u>
Semiquantitative Certification Data:		
I-1	Semiquantitative Certification Testing Analytical Method Conditions	1
I-2	Semiquantitative Certification Testing Statistical Data Summary	2
I-3	DNP - Found Concentrations	3
I-4	DNP - Analysis of Target-Found Concentration Points	4
I-5	RDX - Found Concentrations	6
I-6	RDX - Analysis of Target-Found Concentration Points	7
I-7	TNB - Found Concentrations-HPLC Method	9
I-8	TNB - Analysis of Target-Found Concentration Points- HPLC Method	10
I-9	TNB - Found Concentrations-GC/ECD Method	12
I-10	TNB - Analysis of Target-Found Concentration Points- GC/ECD Method	13
I-11	2,4 - DNT - Found Concentrations-HPLC Method	15
I-12	2,4 - DNT - Analysis of Target-Found Concentration Points-HPLC Method	16
I-13	2,4 - DNT - Found Concentrations-GC/ECD Method	18
I-14	2,4 - DNT - Analysis of Target-Found Concentration Points-GC/ECD Method	19
I-15	TNT - Found Concentrations-HPLC Method	21
I-16	TNT - Analysis of Target-Found Concentration Points- HPLC Method	22

LIST OF TABLES
(Continued)

<u>TABLE NO.</u>		<u>PAGE</u>
I-17	TNT - Found Concentrations-GC/ECD Method	24
I-18	TNT - Analysis of Target Found Concentration Points-GC/ECD Method	25
I-19	Tetryl - Found Concentrations	27
I-20	Tetryl - Analysis of Target-Found Concentration Points	28
I-21	DPA - Found Concentrations	30
I-22	DPA - Analysis of Target-Found Concentration Points	31
I-23	2,6-DNT-Found Concentrations-HPLC Method	33
I-24	2,6-DNT-Analysis of Target-Found Concentration Points-HPLC Method	34
I-25	2,6-DNT-Found Concentrations-GC/ECD Method	36
I-26	2,6-DNT-Analysis of Target-Found Concentration Points-GC/ECD Method	37
I-27	NG-Found Concentrations	39
I-28	NG-Analysis of Target-Found Concentration Points	40
I-29	PETN-Found Concentrations	42
I-30	PETN-Analysis of Target-Found Concentration Points	43

LIST OF TABLES
(Continued)

<u>TABLE NO.</u>		<u>PAGE</u>
	Quantitative Certification Data	
II-1	Quantitative Certification Testing Statistical Data Summary	45
II-2	DNP on Metal-Target vs. Found Concentrations	58
II-3	DNP on Metal-Analysis of Target-Found Concentration Points	59
II-4	DNP on Metal-Inaccuracy and Imprecision Data	61
II-5	DNP on Concrete-Target vs. Found Concentrations	64
II-6	DNP on Concrete-Analysis of Target-Found Concentration Points	65
II-7	DNP on Concrete-Inaccuracy and Imprecision Data	67
II-8	DNP on Concrete (3 days) - Target vs. Found Concentrations	70
II-9	DNP on Concrete (3 days) - Analysis of Target-Found Concentration Points	71
II-10	DNP on Concrete (3 days) - Inaccuracy and Imprecision Data	73
II-11	DNP on Brick - Target vs. Found Concentrations	76
II-12	DNP on Brick - Analysis of Target-Found Concentration Points	77
II-13	DNP on Brick - Inaccuracy and Imprecision Data	79
II-14	DNP on Transite - Target vs. Found Concentrations	82
II-15	DNP on Transite - Analysis of Target-Found Concentration Points	83
II-16	DNP on Transite - Inaccuracy and Imprecision Data	85
II-17	DNP on Transite (3 days) - Target vs. Found Concen- trations	88
II-18	DNP on Transite (3 days) - Analysis of Target-Found Concentration Points	89

LIST OF TABLES
(Continued)

<u>TABLE NO.</u>		<u>PAGE</u>
II-19	DNP on Transite (3 days) - Inaccuracy and Imprecision Data	91
II-20	RDX on Metal - Target vs. Found Concentrations	94
II-21	RDX on Metal - Analysis of Target-Found Concentration Points	95
II-22	RDX on Metal - Inaccuracy and Imprecision Data	97
II-23	RDX on Concrete - Target vs. Found Concentrations	100
II-24	RDX on Concrete - Analysis of Target-Found Concentration Points	101
II-25	RDX on Concrete - Inaccuracy and Imprecision Data	103
II-26	RDX on Brick - Target vs. Found Concentrations	106
II-27	RDX on Brick - Analysis of Target-Found Concentration Points	107
II-28	RDX on Brick - Inaccuracy and Imprecision Data	109
II-29	RDX on Transite - Target vs. Found Concentrations	112
II-30	RDX on Transite - Analysis of Target-Found Concentration Points	113
II-31	RDX on Transite - Inaccuracy and Imprecision Data	115
II-32	RDX on Transite (3 days) - Target vs. Found Concentrations	118
II-33	RDX on Transite (3 days) - Analysis of Target-Found Concentration Points	119
II-34	RDX on Transite (3 days) - Inaccuracy and Imprecision Data	121
II-35	TNB on Metal - Target vs. Found Concentrations	124
II-36	TNB on Metal - Analysis of Target-Found Concentration Points	125

LIST OF TABLES
(Continued)

<u>TABLE NO.</u>		<u>PAGE</u>
II-37	TNB on Metal - Inaccuracy and Imprecision Data	127
II-38	TNB on Concrete - Target vs. Found Concentrations	130
II-39	TNB on Concrete - Analysis of Target-Found Concentration Points	131
II-40	TNB on Concrete - Inaccuracy and Imprecision Data	133
II-41	TNB on Brick - Target vs. Found Concentrations	136
II-42	TNB on Brick - Analysis of Target-Found Concentration Points	137
II-43	TNB on Brick - Inaccuracy and Imprecision Data	139
II-44	TNB on Transite - Target vs. Found Concentrations	142
II-45	TNB on Transite - Analysis of Target-Found Concentration Points	143
II-46	TNB on Transite - Inaccuracy and Imprecision Data	145
II-47	TNB on Transite (3 days) - Target vs. Found Concentrations	148
II-48	TNB on Transite (3 days) - Analysis of Target-Found Concentration Points	149
II-49	TNB on Transite (3 days) - Inaccuracy and Imprecision Data	151
II-50	2,4-DNT on Metal - Target vs. Found Concentrations	154
II-51	2,4-DNT on Metal - Analysis of Target-Found Concentration Points	155
II-52	2,4-DNT on Metal - Inaccuracy and Imprecision Data	157
II-53	2,4-DNT on Concrete - Target vs. Found Concentrations	160
II-54	2,4-DNT on Concrete - Analysis of Target Found Concentration Points	161
II-55	2,4-DNT on Concrete-Inaccuracy and Imprecision Data	163

LIST OF TABLES
(Continued)

<u>TABLE NO.</u>		<u>PAGE</u>
II-56	2,4-DNT on Brick-Target vs. Found Concentrations	166
II-57	2,4-DNT on Brick-Analysis of Target-Found Concentration Points	167
II-58	2,4-DNT on Brick-Inaccuracy and Imprecision Data	169
II-59	2,4-DNT on Transite-Target vs. Found Concentrations	172
II-60	2,4-DNT on Transite-Analysis of Target-Found Concentration Points	173
II-61	2,4-DNT on Transite-Inaccuracy and Imprecision Data	175
II-62	2,4-DNT on Transite (3 days)-Target vs. Found Concentrations	178
II-63	2,4-DNT on Transite (3 days)-Analysis of Target-Found Concentration Points	179
II-64	2,4-DNT on Transite (3 days)-Inaccuracy and Imprecision Data	181
II-65	TNT on Metal - Target vs. Found Concentrations	184
II-66	TNT on Metal - Analysis of Target-Found Concentration Points	185
II-67	TNT on Metal - Inaccuracy and Imprecision Data	187
II-68	TNT on Concrete - Target vs. Found Concentrations	190
II-69	TNT on Concrete - Analysis of Target-Found Concentration Points	191
II-70	TNT on Concrete - Inaccuracy and Imprecision Data	193
II-71	TNT on Brick - Target vs. Found Concentrations	196
II-72	TNT on Brick - Analysis of Target-Found Concentration Points	197
II-73	TNT on Brick - Inaccuracy and Imprecision Data	199
II-74	TNT on Transite - Target vs. Found Concentrations	202

LIST OF TABLES

(Continued)

<u>TABLE NO.</u>		<u>PAGE</u>
II-75	TNT on Transite - Analysis of Target-Found Concentration Points	203
II-76	TNT on Transite - Inaccuracy and Imprecision Data	205
II-77	TNT on Transite (3 days) - Target vs. Found Concentrations	208
II-78	TNT on Transite (3 days) - Analysis of Target-Found Concentration Points	209
II-79	TNT on Transite (3 days) - Inaccuracy and Imprecision Data	211
II-80	Tetryl on Metal - Target vs. Found Concentrations	214
II-81	Tetryl on Metal - Analysis of Target-Found Concentration Points	215
II-82	Tetryl on Metal - Inaccuracy and Imprecision Data	217
II-83	Tetryl on Concrete - Target vs. Found Concentrations	220
II-84	Tetryl on Concrete - Analysis of Target-Found Concentration Points	221
II-85	Tetryl on Concrete - Inaccuracy and Imprecision Data	223
II-86	Tetryl on Brick - Target vs. Found Concentrations	226
II-87	Tetryl on Brick - Analysis of Target-Found Concentration Points	227
II-88	Tetryl on Brick - Inaccuracy and Imprecision Data	229
II-89	Tetryl on Transite - Target vs. Found Concentrations	232
II-90	Tetryl on Transite - Analysis of Target-Found Concentration Points	233
II-91	Tetryl on Transite - Inaccuracy and Imprecision Data	235
II-92	DPA on Metal - Target vs. Found Concentrations	238
II-93	DPA on Metal - Analysis of Target-Found Concentration Points	239

LIST OF TABLES
(Continued)

<u>TABLE NO.</u>		<u>PAGE</u>
II-94	DPA on Metal - Inaccuracy and Imprecision Data	241
II-95	DPA on Concrete - Target vs. Found Concentrations	244
II-96	DPA on Concrete - Analysis of Target-Found Concentration Points	245
II-97	DPA on Concrete - Inaccuracy and Imprecision Data	247
II-98	DPA on Concrete (3 days) - Target vs. Found Concentrations	250
II-99	DPA on Concrete (3 days) - Analysis of Target-Found Concentration Points	251
II-100	DPA on Concrete (3 days) - Inaccuracy and Imprecision Data	253
II-101	DPA on Brick - Target vs. Found Concentrations	256
II-102	DPA on Brick - Analysis of Target-Found Concentration Points	257
II-103	DPA on Brick - Inaccuracy and Imprecision Data	259
II-104	DPA on Transite - Target vs. Found Concentrations	262
II-105	DPA on Transite - Analysis of Target-Found Concentration Points	263
II-106	DPA on Transite - Inaccuracy and Imprecision Data	265
II-107	DPA on Transite (3 days) - Target vs. Found Concentrations	268
II-108	DPA on Transite (3 days) - Analysis of Target-Found Concentration Points	269
II-109	DPA on Transite (3 days) - Inaccuracy and Imprecision Data	271
II-110	2,6-DNT on Metal - Target vs. Found Concentrations	274
II-111	2,6-DNT on Metal - Analysis of Target-Found Concentration Points	275

LIST OF TABLES
(Continued)

<u>TABLE NO.</u>		<u>PAGE</u>
II-112	2,6-DNT on Metal - Inaccuracy and Imprecision Data	277
II-113	2,6-DNT on Concrete - Target vs. Found Concentrations	280
II-114	2,6-DNT on Concrete - Analysis of Target-Found Concentration Points	281
II-115	2,6-DNT on Concrete - Inaccuracy and Imprecision Data	283
II-116	2,6-DNT on Brick - Target vs. Found Concentrations	286
II-117	2,6-DNT on Brick - Analysis of Target-Found Concentration Points	287
II-118	2,6-DNT on Brick - Inaccuracy and Imprecision Data	289
II-119	2,6-DNT on Transite - Target vs. Found Concentrations	292
II-120	2,6-DNT on Transite - Analysis of Target-Found Concentration Points	293
II-121	2,6-DNT on Transite - Inaccuracy and Imprecision Data	295
II-122	NG on Metal - Target vs. Found Concentrations	298
II-123	NG on Metal - Analysis of Target-Found Concentration Points	299
II-124	NG on Metal - Inaccuracy and Imprecision Data	301
II-125	NG on Concrete Target vs. Found Concentrations	304
II-126	NG on Concrete - Analysis of Target-Found Concentration Points	305
II-127	NG on Concrete - Inaccuracy and Imprecision Data	307
II-128	NG on Brick - Target vs. Found Concentrations	310
II-129	NG on Brick - Analysis of Target-Found Concentration Points	311

LIST OF TABLES
(Continued)

<u>TABLE NO.</u>		<u>PAGE</u>
II-130	NG on Brick - Inaccuracy and Imprecision Data	313
II-131	NG on Transite - Target vs. Found Concentrations	316
II-132	NG on Transite - Analysis of Target-Found Concentration Points	317
II-133	NG on Transite - Inaccuracy and Imprecision Data	319
II-134	PETN on Metal - Target vs. Found Concentrations	322
II-135	PETN on Metal - Analysis of Target-Found Concentration Points	324
II-136	PETN on Concrete - Target vs. Found Concentrations	326
II-137	PETN on Concrete - Analysis of Target-Found Concentration Points	328
II-138	PETN on Brick - Target vs. Found Concentrations	33-
II-139	PETN on Brick - Analysis of Target-Found Concentration Points	332
II-140	PETN on Transite - Target vs. Found Concentrations	334
II-141	PETN on Transite - Analysis of Target-Found Concentration Points	336

LIST OF FIGURES

<u>FIGURE NO.</u>		<u>PAGE</u>
Semiquantitative Certification Data:		
I-1	DNP - Graph of Target-Found Concentration Points	5
I-2	RDX - Graph of Target-Found Concentration Points	8
I-3	TNB - Graph of Target-Found Concentration Points- HPLC Method	11
I-4	TNB - Graph of Target-Found Concentration Points- GC/ECD Method	14
I-5	2,4-DNT-Graph of Target-Found Concentration Points-HPLC Method	17
I-6	2,4-DNT-Graph of Target-Found Concentration Points-GC/ECD Method	20
I-7	TNT - Graph of Target-Found Concentration Points-HPLC Method	25
I-8	TNT - Graph of Target-Found Concentration Points-GC/ECD Method	26
I-9	Tetryl - Graph of Target-Found Concentration Points	29
I-10	DPA - Graph of Target-Found Concentration Points	32
I-11	2,6-DNT-Graph of Target-Found Concentration Points- HPLC Method	35
I-12	2,6-DNT-Graph of Target-Found Concentration Points- GC/ECD Method	38
I-13	NG - Graph of Target-Found Concentration Points	41
I-14	PETN - Graph of Target-Found Concentration Points	44

LIST OF FIGURES

<u>FIGURE NO.</u>		<u>PAGE</u>
	Quantitative Certification Data	
II-1	DNP on Metal - Graph of Target-Found Concentration Points	60
II-2	DNP on Metal - Graph of Inaccuracy	62
II-3	DNP on Metal - Graph of Imprecision	63
II-4	DNP on Concrete - Graph of Target-Found Concentration Points	66
II-5	DNP on Concrete - Graph of Inaccuracy	68
II-6	DNP on Concrete - Graph of Imprecision	69
II-7	DNP on Concrete (3 days) - Graph of Target-Found Concentration Points	72
II-8	DNP on Concrete (3 days) - Graph of Inaccuracy	74
II-9	DNP on Concrete (3 days) - Graph of Imprecision	75
II-10	DNP on Brick - Graph of Target-Found Concentration Points	78
II-11	DNP on Brick - Graph of Inaccuracy	80
II-12	DNP on Brick - Graph of Imprecision	81
II-13	DNP on Transite - Graph of Target-Found Concentration Points	84
II-14	DNP on Transite - Graph of Inaccuracy	86
II-15	DNP on Transite - Graph of Imprecision	87
II-16	DNP on Transite (3 days) - Graph of Target-Found Concentration Points	90
II-17	DNP on Transite (3 days) - Graph of Inaccuracy	92
II-18	DNP on Transite (3 days) - Graph of Imprecision	93
II-19	RDX on Metal - Graph of Target-Found Concentration Points	96

LIST OF FIGURES
(Continued)

<u>FIGURE NO.</u>		<u>PAGE</u>
II-20	RDX on Metal - Graph of Inaccuracy	98
II-21	RDX on Metal - Graph of Imprecision	99
II-22	RDX on Concrete - Graph of Target-Found Concentration Points	102
II-23	RDX on Concrete - Graph of Inaccuracy	104
II-24	RDX on Concrete - Graph of Imprecision	105
II-25	RDX on Brick - Graph of Target-Found Concentration Points	108
II-26	RDX on Brick - Graph of Inaccuracy	110
II-27	RDX on Brick - Graph of Imprecision	111
II-28	RDX on Transite - Graph of Target-Found Concentration Points	114
II-29	RDX on Transite - Graph of Inaccuracy	116
II-30	RDX on Transite - Graph of Imprecision	117
II-31	RDX on Transite (3 days) - Graph of Target-Found Concentration Points	120
II-32	RDX on Transite (3 days) - Graph of Inaccuracy	122
II-33	RDX on Transite (3 days) - Graph of Imprecision	123
II-34	TNB on Metal - Graph of Target-Found Concentration Points	126
II-35	TNB on Metal - Graph of Inaccuracy	128
II-36	TNB on Metal - Graph of Imprecision	129
II-37	TNB on Concrete - Graph of Target-Found Concentration Points	132
II-38	TNB on Concrete - Graph of Inaccuracy	134
II-39	TNB on Concrete - Graph of Imprecision	135

LIST OF FIGURES
(Continued)

<u>FIGURE NO.</u>		<u>PAGE</u>
II-40	TNB on Brick - Graph of Target-Found Concentration Points	138
II-41	TNB on Brick - Graph of Inaccuracy	140
II-42	TNB on Brick - Graph of Imprecision	141
II-43	TNB on Transite - Graph of Target-Found Concentration Points	144
II-44	TNE on Transite - Graph of Inaccuracy	146
II-45	TNB on Transite - Graph of Imprecision	147
II-46	TNB on Transite (3 days) - Graph of Target-Found Concentration Points	150
II-47	TNB on Transite (3 days) - Graph of Inaccuracy	152
II-48	TNB on Transite (3 days) - Graph of Imprecision	153
II-49	2,4-DNT on Metal - Graph of Target-Found Concentration Points	156
II-50	2,4-DNT on Metal - Graph of Inaccuracy	158
II-51	2,4-DNT on Metal - Graph of Imprecision	159
II-52	2,4-DNT on Concrete - Graph of Target-Concentration Points	162
II-53	2,4-DNT on Concrete - Graph of Inaccuracy	164
II-54	2,4-DNT on Concrete - Graph of Imprecision	165
II-55	2,4-DNT on Brick - Graph of Target-Found Concentration Points	168
II-56	2,4-DNT on Brick - Graph of Inaccuracy	170
II-57	2,4-DNT on Brick - Graph of Imprecision	171
II-58	2,4-DNT on Transite - Graph of Target-Found Concentration Points	174

LIST OF FIGURES
(Continued)

<u>FIGURE NO.</u>		<u>PAGE</u>
II-59	2,4-DNT on Transite - Graph of Inaccuracy	176
II-60	2,4-DNT on Transite - Graph of Imprecision	177
II-61	2,4-DNT on Transite (3 days) - Graph of Target- Found Concentration Points	180
II-62	2,4-DNT on Transite (3 days) - Graph of Inaccuracy	182
II-63	2,4-DNT on Transite (3 days) - Graph of Imprecision	183
II-64	TNT on Metal - Graph of Target-Found Concentration Points	186
II-65	TNT on Metal - Graph of Inaccuracy	188
II-66	TNT on Metal - Graph of Imprecision	189
II-67	TNT on Concrete - Graph of Target-Found Concen- tration Points	192
II-68	TNT on Concrete - Graph of Inaccuracy	194
II-69	TNT on Concrete - Graph of Imprecision	195
II-70	TNT on Brick - Graph of Target-Found Concen- tration Points	198
II-71	TNT on Brick - Graph of Inaccuracy	200
II-72	TNT on Brick - Graph of Imprecision	201
II-73	TNT on Transite - Graph of Target-Found Concen- tration Points	204
II-74	TNT on Transite - Graph of Inaccuracy	206
II-75	TNT on Transite - Graph of Imprecision	207
II-76	TNT on Transite (3 days) - Graph of Target-Found Concentration Points	210
II-77	TNT on Transite (3 days) - Graph of Inaccuracy	212
II-78	TNT on Transite (3 days) - Graph of Imprecision	213

LIST OF FIGURES
(Continued)

<u>FIGURE NO.</u>		<u>PAGE</u>
II-79	Tetryl on Metal - Graph of Target-Found Concentration Points	216
II-80	Tetryl on Metal - Graph of Inaccuracy	218
II-81	Tetryl on Metal - Graph of Imprecision	219
II-82	Tetryl on Concrete - Graph of Target-Found Concentration Points	222
II-83	Tetryl on Concrete - Graph of Inaccuracy	224
II-84	Tetryl on Concrete - Graph of Imprecision	225
II-85	Tetryl on Brick - Graph of Target-Found Concentration Points	228
II-86	Tetryl on Brick - Graph of Inaccuracy	230
II-87	Tetryl on Brick - Graph of Imprecision	231
II-88	Tetryl on Transite - Graph of Target-Found Concentration Points	234
II-89	Tetryl on Transite - Graph of Inaccuracy	236
II-90	Tetryl on Transite - Graph of Imprecision	237
II-91	DPA on Metal - Graph of Target-Found Concentration Points	240
II-92	DPA on Metal - Graph of Inaccuracy	242
II-93	DPA on Metal - Graph of Imprecision	243
II-94	DPA on Concrete - Graph of Target-Found Concentration Points	246
II-95	DPA on Concrete - Graph of Inaccuracy	248
II-96	DPA on Concrete - Graph of Imprecision	249
II-97	DPA on Concrete (3 days) - Graph of Target-Found Concentration Points	252

LIST OF FIGURES
(Continued)

<u>FIGURE NO.</u>		<u>PAGE</u>
II-98	DPA on Concrete (3 days) - Graph of Inaccuracy	254
II-99	DPA on Concrete (3 days) - Graph of Imprecision	255
II-100	DPA on Brick - Graph of Target-Found Concentration Points	258
II-101	DPA on Brick - Graph of Inaccuracy	260
II-102	DPA on Brick - Graph of Imprecision	261
II-103	DPA on Transite - Graph of Target-Found Concentration Points	264
II-104	DPA on Transite - Graph of Inaccuracy	266
II-105	DPA on Transite - Graph of Imprecision	267
II-106	DPA on Transite (3 days) - Graph of Target-Found Concentration Points	270
II-107	DPA on Transite (3 days) - Graph of Inaccuracy	272
II-108	DPA on Transite (3 days) - Graph of Imprecision	273
II-109	2,6-DNT on Metal - Graph of Target-Found Concentration Points	276
II-110	2,6-DNT on Metal - Graph of Inaccuracy	278
II-111	2,6-DNT on Metal - Graph of Imprecision	279
II-112	2,6-DNT on Concrete - Graph of Target-Found Concentration Points	282
II-113	2,6-DNT on Concrete - Graph of Inaccuracy	284
II-114	2,6-DNT on Concrete - Graph of Imprecision	285
II-115	2,6-DNT on Brick - Graph of Target-Found Concentration Points	288
II-116	2,6-DNT on Brick - Graph of Inaccuracy	290
II-117	2,6-DNT on Brick - Graph of Imprecision	291

LIST OF FIGURES
(Continued)

<u>FIGURE NO.</u>		<u>PAGE</u>
II-118	2,6-DNT on Transite - Graph of Target-Found Concentration Points	294
II-119	2,6-DNT on Transite - Graph of Inaccuracy	296
II-120	2,6-DNT on Transite - Graph of Imprecision	297
II-121	NG on Metal - Graph of Target-Found Concentration Points	300
II-122	NG on Metal - Graph of Inaccuracy	302
II-123	NG on Metal - Graph of Imprecision	303
II-124	NG on Concrete - Graph of Target-Found Concentration Points	306
II-125	NG on Concrete - Graph of Inaccuracy	308
II-126	NG on Concrete - Graph of Imprecision	309
II-127	NG on Brick - Graph of Target-Found Concentration Points	312
II-128	NG on Brick - Graph of Inaccuracy	314
II-129	NG on Brick - Graph of Imprecision	315
II-130	NG on Transite - Graph of Target-Found Concentration Points	318
II-131	NG on Transite - Graph of Inaccuracy	320
II-132	NG on Transite - Graph of Imprecision	321
II-133	PETN on Metal - Graph of Target-Found Concentration Points	325
II-134	PETN on Concrete - Graph of Target-Found Concentration Points	329
II-135	PETN on Brick - Graph of Target-Found Concentration Points	333
II-136	PETN on Transite - Graph of Target-Found Concentration Points	337

Table I-1. Semiquantitative Certification Testing!
Analytical Method Conditions

Analyte	Instrument	Column	Precolumn	Program	Retention Time	Solvent System/ Carrier Gas	Flow Rate
PEIN	HPIC	Spherisorb ODS, 5µ, 250 x 4.6 mm	Pellicular LC-18, 40µ, 50 x 4.6 mm	Isocratic	6.2 min	65/35 CH ₃ CH ₃ /H ₂ O	1 mL/min
2,6-DNT	HPIC	Spherisorb ODS, 5µ, 250 x 4.6 mm	Pellicular LC-18, 40µ, 50 x 4.6 mm	Isocratic	24.6 min	35/65 CH ₃ CN/ .005M t-butyl	1 mL/min
	HPIC	Spherisorb ODS, 5µ, 250 x 4.6 mm	Pellicular LC-18 40µ, 50 x 4.6 mm	Isocratic	30.5 min	ammonium hydroxide	1 mL/min
	GC	3% OV-225 on 100/120 Gas ChromQ 1/8" x 2 mm ID x 6' glass column	None	100C for 6 min 15C/min to 165C _i	19.90 min	5% Methane/argon	30 mL/min
1,3,5-TNB	GC	3% OV-225 on 100/120 Gas ChromQ 1/8" x 2 mm ID x 6' glass column	None	hold 8 min 15C/min to 200 C _i	16.20 min	5% Methane/argon	30 mL/min
	GC	3% OV-225 on 100/120 Gas ChromQ 1/8" x 2 mm ID x 6' glass column	None	hold 6 min	30.95 min	5% Methane/argon	30 mL/min
	GC	3% OV-225 on 100/120 Gas ChromQ 1/8" x 2 mm ID x 6' glass column	None	"	28.14 min	5% Methane/argon	30 mL/min
DPA	HPIC	U Bondpack C18 4mm x 30 cm	None	Isocratic	2.30 sec	90/10 methanol/H ₂ O	1 mL/min
Tetryl	HPIC	U Bondpack C18 4mm x 30 cm	None	Isocratic	416 sec	60/40 methanol/H ₂ O	1 mL/min
DNT	HPIC	Spherisorb ODS, 5µ, 250 x 4.6 mm	Pellicular LC-18, 40µ, 50 x 4.6 mm	Initial: 30% CH ₃ CN Final: 50% CH ₃ CN Time: 35 min Gradient: linear	12.5 min	0.08M acetic acid adjusted to pH 3.1 with ammonium hydroxide/CH ₃ CN	1.0 mL/min
RDX	HPIC	Spherisorb ODS, 5µ, 250 x 4.6 mm	Pellicular LC-18, 40µ, 50 x 4.6 mm	Initial: 30% CH ₃ CN Final: 50% CH ₃ CN Time: 35 min Gradient: linear	13.5 min	0.08M acetic acid adjusted to pH 3.1 with ammonium hydroxide/CH ₃ CN	1.0 mL/min
TNB	HPIC	Spherisorb ODS, 5µ, 250 x 4.6 mm	Pellicular LC-18, 40µ, 50 x 4.6 mm	Initial: 30% CH ₃ CN Final: 50% CH ₃ CN Time: 35 min Gradient: linear	17.2 min	0.08M acetic acid adjusted to pH 3.1 with ammonium hydroxide/CH ₃ CN	1.0 mL/min
2,4-DNT	HPIC	Spherisorb ODS, 5µ, 250 x 4.6 mm	Pellicular LC-18, 40µ, 50 x 4.6 mm	Initial: 30% CH ₃ CN Final: 50% CH ₃ CN Time: 35 min Gradient: linear	24.3 min	0.08M acetic acid adjusted to pH 3.1 with ammonium hydroxide/CH ₃ CN	1.0 mL/min
TNT	HPIC	Spherisorb ODS, 5µ, 250 x 4.6 mm	Pellicular LC-18, 40µ, 50 x 4.6 mm	Initial: 30% CH ₃ CN Final: 50% CH ₃ CN Time: 35 min Gradient: linear	25.6 min	0.08M acetic acid adjusted to pH 3.1 with ammonium hydroxide/CH ₃ CN	1.0 mL/min
Tetryl	HPIC	Spherisorb ODS, 5µ, 250 x 4.6 mm	Pellicular LC-18, 40µ, 50 x 4.6 mm	Initial: 30% CH ₃ CN Final: 50% CH ₃ CN Time: 35 min Gradient: linear	26.7 min	0.08M acetic acid adjusted to pH 3.1 with ammonium hydroxide/CH ₃ CN	1.0 mL/min
DPA	HPIC	Spherisorb ODS, 5µ, 250 x 4.6 mm	Pellicular LC-18, 40µ, 50 x 4.6 mm	Initial: 30% CH ₃ CN Final: 50% CH ₃ CN Time: 35 min Gradient: linear	39.2 min	0.08M acetic acid adjusted to pH 3.1 with ammonium hydroxide/CH ₃ CN	1.0 mL/min

Table I-2. Semiquantitative Certification Testing
Statistical Data Summary

Analyte	Detector	Attn.	Chart Speed	Inj. Volume	USAT/AMA Method	Det. Limit	Corr. Coeff.	Slope	Int.	MR Reference
PEIN	UV at 230 nm	0.01 AUFS	0.1 in/min	70 uL	none	1.77 ug/mL	0.996	1.046	0.450	9
{ 2,6-DNT NI	UV at 230 nm	0.01 AUFS	0.1 in/min	70 uL	MRI Method	0.26 ug/mL	0.998	1.061	9.888	8
	UV at 230 nm	0.01 AUFS	0.1 in/min	70 uL	MRI Method	4.54 ug/mL	0.999	0.989	2119.5	8
	ECD @300C	1x10 ⁻¹¹ g	0.5 cm/min	1 uL	246 TNT-MA-02	0.11 ug/mL	0.967	0.743	0.034	5
{ 2,6-DNT 1,3,5-TNB 2,4,6-TNT	ECD @300C	1x10 ⁻¹¹ g	0.5 cm/min	1 uL	246 TNT-MA-02	0.09 ug/mL	0.989	1.069	-0.012	5
	ECD @300C	1x10 ⁻¹¹ g	0.5 cm/min	1 uL	246 TNT-MA-02	0.12 ug/mL	0.999	1.058	0.007	5
	ECD @300C	1x10 ⁻¹¹ g	0.5 cm/min	1 uL	246 TNT-MA-07	0.12 ug/mL	0.979	1.054	-0.009	5
DPA	UV at 254 nm	0.1 AUFS	0.5 cm/min	200 uL	DPA-MA-01	12.25 ug/mL	0.997	0.796	-0.634	3
Tetryl	UV at 254 nm	0.1 AUFS	0.5 cm/min	200 uL	TETRYL-MA-02	0.08 ug/mL	0.995	0.862	-0.004	5
UMP	UV at 254 nm	0.01 AUFS	0.1 in/min	100 uL	MRI Method	30 ng/mL	0.998	0.971	5.397	10
BOX	UV at 254 nm	0.01 AUFS	0.1 in/min	100 uL	MIR Method	70 ng/mL	0.998	1.009	9.698	10
TNB	UV at 254 nm	0.01 AUFS	0.1 in/min	100 uL	MRI Method	25 ng/mL	0.999	0.994	0.674	10
{ 2,4-DNT TNT	UV at 254 nm	0.01 AUFS	0.1 in/min	100 uL	MRI Method	24 ng/mL	0.999	0.990	1.585	10
	UV at 254 nm	0.01 AUFS	0.1 in/min	100 uL	MRI Method	24 ng/mL	0.999	1.012	0.269	10
Tetryl	UV at 254 nm	0.01 AUFS	0.1 in/min	100 uL	MRI Method	25 ng/mL	0.999	1.050	-5.410	10
DPA	UV at 254 nm	0.01 AUFS	0.1 in/min	100 uL	MRI Method	51 ng/mL	0.999	0.993	7.745	10

Table I-3. DNP - Found Concentrations

ARMY CERTIFICATION	
DNP	FOUND CONC
(1)	26 24 (2) 26 28
(3)	26 29 (4) 52 52
(5)	52 52 (6) 52 53
(7)	104 104 (8) 104 126
(9)	104 111 (10) 260 254
(11)	260 261 (12) 260 264
(13)	520 495 (14) 520 533
(15)	520 498

Table I-4. DNP - Analysis of Target-Found Concentration Points

ARMY CERTIFICATION
DNP
ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 192.4 SD 189.231830605

FOUND CONC
MEAN= 192.266666667 SD= 184.077804461

NO. RUNS : TOTAL X-Y ALL RUNS 15 NO. CONCENTR 15
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 5.39670981661
SLOPE= 0.971257571986
USE FOR ACCURACY
R= 0.998451979989
MEAN SQR DEV OF POINTS FROM REGRESSION= 112.090610730
ST ERROR EST= 10.6249993288
USE FOR PRECISION
T FOR CONFIDENCE BAND
D.F.= 13
TWO TAIL P LEVEL IS .1
t= 1.77093170942
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U U
REPLICATES= 2
y(c)= 19.6057408586
x(d)= 29.1259198183

ARMY CERTIFICATION
DNP
FOUND CONC

533.00

$y(c)$

$x(d)$

0.00

VERTICAL AXIS TIC INTERVAL = 53.3
TARGET CONC

520

416

312

208

104

Figure I-1. DNP - Graph of Target-Found Concentration Points

Table I-5. RDX - Found Concentrations

ARMY CERTIFICATION	
RDX	FOUND CONC
(1)	50 44 (2) 50 44
(3)	50 60 (4) 100 95
(5)	100 107 (6) 100 110
(7)	200 235 (8) 200 193
(9)	200 271 (10) 500 500
(11)	500 516 (12) 500 521
(13)	1000 1000 (14) 1000 1066
(15)	1000 984

Table I-6. RDX - Analysis of Target-Found Concentration Points

ARMY CERTIFICATION
RDX
ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 378 SD= 363.907366349

FOUND CONC
MEAN= 383.066666667 SD= 368.11963635

NO. RUNS 1 TOTAL X-Y ALL RUNS 15 NO. CONCENTR 15
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 9.69795037756
SLOPE= 1.00910463862
USE FOR ACCURACY
R= 0.997557791953
MEAN SQR DEV OF POINTS FROM REGRESSION= 711.942079499
ST ERROR EST= 26.6822427749
USE FOR PRECISION
T FOR CONFIDENCE BAND
D.F.= 13
TWO TAIL P LEVEL IS .1
t= 1.77093170942
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U U
REPLICATES= 2
u(c)= 45.3006048026
x(d)= 70.3255803599

ARMY CERTIFICATION

RDX

FOUND CONC

1066.00




Figure I-2. RDX - Graph of Target-Found Concentration Points

Table I-7. TNB - Found Concentrations - HPLC Method

ARMY CERTIFICATION

TNB

FOUND CONC

<1> 25 25 <2> 25 26

<3> 25 28 <4> 50 49

<5> 50 50 <6> 50 52

<7> 100 96 <8> 100 99

<9> 100 104 <10> 250 243

<11> 250 253 <12> 250 250

<13> 500 486 <14> 500 518

<15> 500 491

Table I-8. TNE - Analysis of Target-Found
Concentration Points - HPLC Method

ARMY CERTIFICATION
TNE
ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 185 SD= 181.953683274

FOUND CONC
MEAN= 184.666666667 SD= 181.09889432

N0. RUNS 1 TOTAL X-Y ALL RUNS 15 N0. CONCENTR 15
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.674487594391
SLOPE= 0.99455231931
USE FOR ACCURACY
R= 0.999246618194
MEAN SQR DEV OF POINTS FROM REGRESSION= 53.1983838461
ST ERROR EST= 7.2937166826
USE FOR PRECISION
T FOR CONFIDENCE BAND
D.F.= 13
TWO TAIL P LEVEL IS .1
t= 1.77893178942
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U U
REPLICATES= 2
y(c)= 10.4285310439
x(d)= 19.5311771834

ARMY CERTIFICATION

TNB

FOUND CONC

518.00




Figure 1-3. TNB - Graph of Target-Found Concentration Points - HPLC Method

Table I-9. TNB-Found Concentrations - GC/ECD Method

ARMY CERTIFICATION GC/ECD METHOD

1,3,5 TNB

FOUND CONC

(1)	0.12	0.09	(2)	0.12	0.1
(3)	0.12	0.09	(4)	0.12	0.17
(5)	0.23	0.22	(6)	0.23	0.31
(7)	0.23	0.25	(8)	0.46	0.53
(9)	0.46	0.55	(10)	0.46	0.48
(11)	0.46	0.49	(12)	1.16	1.23
(13)	1.16	1.27	(14)	1.16	1.26
(15)	2.31	2.42	(16)	2.31	2.44
(17)	2.31	2.46	(18)	3	0
(19)	0	0	(20)	0	0
(21)	0	0			

Table I-10. TNB - Analysis of Target-Found Concentration Points - GC/ECD Method

ARMY CERTIFICATION GC/ECD METHOD

1.3.5 TNB

ANALYSIS OF 21 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 0.639047619048 SD= 0.791573779812

FOUND CONC

MEAN= 0.68380952381 SD= 0.838793634874

NB. RUNS 1 TOTAL X-Y ALL RUNS 21 NB. CONCENTR 21
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.60709275380941

SLOPE= 1.05894575030

R= 0.99332439997

MEAN SQR DEV OF POINTS FROM REGRESSION (ST ERROR EST)= 3.884665285E-4
COMPUTE T

D.F. = 19

ENTER 2 TAIL P LEVEL (USUALLY .1. EACH CONFID BAND IS

.05 SO TOTAL P= .1)

.1

t= 1.729129474

REPLICATES ON UNKNOWN SAMPLE 3

Y(CC)= 0.0420524389347

X(CC)= 0.0657723484472

ARMY CERTIFICATION GC/ECD METHOD
 1,3,5 TNB
 FOUND CONC




Figure I-4. TNB - Graph of Target-Found Concentration Points - GC/ECD Method

Table I-11. 2,4-DNT - Found Concentrations - HPLC Method

ARMY CERTIFICATION

24DNT

FOUND CONC

(1) 24 26 (2) 24 26

(3) 24 27 (4) 48 50

(5) 48 48 (6) 48 50

(7) 96 91 (8) 96 94

(9) 96 100 (10) 240 238

(11) 240 241 (12) 240 240

(13) 480 471 (14) 480 492

(15) 480 468

Table I-12. 2,4-DNT - Analysis of Target-Found
Concentration Points - HPLC Method

ARMY CERTIFICATION
24DNT
ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 177.6 SD= 174.675535943

FOUND CONC
MEAN= 177.466666667 SD= 173.068882516

NO. RUNS 1 TOTAL X-Y ALL RUNS 15 NO. CONCENTR 15
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 1.58495145631
SLOPE= 0.990324973031
USE FOR ACCURACY
R= 0.999518474424
MEAN SQR DEV OF POINTS FROM REGRESSION= 31.037567837
ST ERROR EST= 5.37293170935
USE FOR PRECISION
T FOR CONFIDENCE BAND
D.F. = 13
TWO TAIL P LEVEL IS .1
t= 1.77093170942
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U U
REPLICATES= 2
y(c)= 9.03775232583
x(d)= 15.0116760434

ARMY CERTIFICATION

24DNT
FOUND CONC

492.00

0

y(c)

x(d)

0
0.00

TARGET CONC

VERTICAL AXIS TIC INTERVAL = 49.2

480

384

288

192

96

Figure I-5. 2,4-DNT - Graph of Target-Found Concentration
Points - HPLC Method

Table I-13. 2,4-DNT-Found Concentrations - GC/ECD Method

ARMY CERTIFICATION GC/ECD METHOD

2,4 DNT

FOUND CONC

(1)	0.03	0.05	(2)	0.03	0.04
(3)	0.03	0.06	(4)	0.06	0.06
(5)	0.06	0.08	(6)	0.06	0.08
(7)	0.11	0.1	(8)	0.11	0.11
(9)	0.11	0.14	(10)	0.11	0.21
(11)	0.28	0.2	(12)	0.28	0.27
(13)	0.28	0.31	(14)	0.28	0.31
(15)	0.56	0.43	(16)	0.56	0.42
(17)	0.56	0.42	(18)	0.03	0.07
(19)	0	0	(20)	0	0
(21)	0	0.01	(22)	0	0

Table I-14. 2,4-DNT - Analysis of Target-Found Concentration Points - GC/ECD Method

ARMY CERTIFICATION GC/ECD METHOD
 2,4 DNT
 ANALYSIS OF 22 TARGET CONC-FOUND CONC POINTS

TARGET CONC
 MEAN= 0.160909090909 SD= 0.188677328259

FOUND CONC
 MEAN= 0.153181818182 SD= 0.144992162798

NO. RUNS 1 TOTAL X-Y ALL RUNS 22 NO. CONCENTR 22
 MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.0336084831092
 SLOPE= 0.7431112192
 R= 0.367805515185
 MEAN SQ. DEV OF POINTS FROM REGRESSION (ST ERROR EST)= 0.0014326811139
 COMPUTE 1

D.F. = 20
 ENTER 2 TAIL P LEVEL (USUALLY .1, EACH CONFID BAND IS
 .05 SO TOTAL P= .1)
 .1
 t= 1.72471716621
 REPLICATES ON UNKNOWN SAMPLE 3
 y(c)= 0.0755821824279
 x(d)= 0.110789286561

Figure 1 is a graph showing the relationship between Target Concentration (X-axis) and TIC Interval (Y-axis). The X-axis ranges from 0.00 to 0.56, and the Y-axis ranges from 0.00 to 0.13. Two parallel lines are plotted, with data points marked by circles. A vertical line is drawn at X=0.112, intersecting the lower line at Y=0.043. A horizontal line is drawn at Y=0.043, intersecting the upper line at X=0.112. The vertical axis is labeled 'VERTICAL AXIS TIC INTERVAL = 0.043'.

Table I-15. TNT - Found Concentrations - HPLC Method

ARMY CERTIFICATION	
TNT	
FOUND CONC	
(1) 24 25 (2) 24 26	
(3) 24 26 (4) 48 48	
(5) 48 49 (6) 48 50	
(7) 96 89 (8) 96 98	
(9) 96 102 (10) 244 242	
(11) 244 248 (12) 244 250	
(13) 488 490 (14) 488 507	
(15) 488 487	

Table I-16. TNT - Analysis of Target-Found Concentration
Points - HPLC Method

ARMY CERTIFICATION

TNT

ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 180 SD= 177.950233014

FOUND CONC
MEAN= 182.46666667 SD= 180.195475869

NO. RUNS 1 TOTAL X-Y ALL RUNS 15 NO. CONCENTR 15
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.269286848564
SLOPE= 1.01220766566
USE FOR ACCURACY
R= 0.999595522042
MEAN SQR DEV OF POINTS FROM REGRESSION= 28.2819574468
ST ERROR EST= 5.31807836035

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 13

TWO TAIL P LEVEL IS .1

t= 1.77093170942

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLET C/U U

REPLICATES= 2

y(c)= 7.3789069858

x(d)= 14.0139688505




Figure I-7. TNT - Graph of Target-Found Concentration Points - HPLC Method

Table I-17. TNT - Found Concentrations - GC/ECD Method

ARMY CERTIFICATION GC/ECD METHOD:

2,4,6 TNT

FOUND CONC

(1) 0.04 0.02 (2) 0.04 0.03

(3) 0.04 0.1 (4) 0.04 0.00

(5) 0.00 0.12 (5) 0.00 0

(7) 0.00 0.06 (8) 0.00 0.15

(9) 0.15 0.11 (10) 0.15 0.14

(11) 0.15 0.00 (12) 0.15 0.00

(13) 0.37 0.51 (14) 0.37 0.32

(15) 0.37 0.34 (16) 0.37 0.34

(17) 0.75 0.92 (18) 0.75 0.77

(19) 0.75 0.73 (20) 0.75 0.74

(21) 0 0 (22) 0 0

(23) 0 0 (24) 0 0

Table I-18. TNT - Analysis of Target Found Concentration
Points - GC/ECD Method

ARMY CERTIFICATION GC/ECD METHOD

2,4,6 TNT

ANALYSIS OF 24 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 0.23166666667 SD= 0.266338566272

FOUND CONC

MEAN= 0.235 SD= 0.28672858121

NB. RUNS 1 TOTAL X-Y ALL RUNS 24 NO. CONCENTR 24
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.08911473269162

SLOPE= 1.053732685

R= 0.978826309535

MEAN SQR DEV OF POINTS FROM REGRESSION (ST ERROR EST)= 0.00360183498674
COMPUTE T

D.F.= 22

ENTER 2 TAIL P LEVEL (USUALLY .1) CONFID BAND IS

.05 SO TOTAL P= .1

.1

t= 1.71713909197

REPLICATES ON UNKNOWN SAMPLE 3

y(c)= 0.0566954873641

x(d)= 0.122913743297

ARMY CERTIFICATION GC/ECD METHOD
2,4,6 TNT
FOUND CONC




Figure 1-8. Graph of Target-Found Concentration
Points - GC/ECD Method

Table I-19. Tetryl-Found Concentrations

ARMY CERTIFICATION	
TETRYL	
FOUND CONC	
RUN 2	
(14) 0 0 (15) 0 0	
(16) 0 0 (17) 0 0	
(18) 0.05 0.042 (19) 0.05 0.039	
(20) 0.05 0.039 (21) 0.05 0.036	
(22) 0.1 0.086 (23) 0.1 0.086	
(24) 0.1 0.083 (25) 0.1 0.083	
(26) 0.2 0.172 (27) 0.2 0.179	
(28) 0.2 0.166 (29) 0.2 0.256	
(30) 0.5 0.367 (31) 0.5 0.383	
(32) 0.5 0.364 (33) 0.5 0.37	
(34) 1 0.874 (35) 1 0.888	
(36) 1 0.872 (37) 1 0.885	

Table I-20. Tetryl - Analysis of Target-Found Concentration Points

ARMY CERTIFICATION

TETRYL

ANALYSIS OF 24 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 0.308333333333 SD= 0.356817657326

FOUND CONC

MEAN= 0.26125 SD= 0.309147108643

N0. RUNS 1 TOTAL X-Y ALL RUNS 24 N0. CONCENTR 24
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.00450455321571

SLOPE= 0.861906659078

R= 0.994812845018

MEAN SQR DEV OF POINTS FROM REGRESSION (ST ERROR EST)= 0.00103307217364
COMPUTE 1

D.F. = 22

ENTER 2 TAIL P LEVEL (USUALLY .1, EACH CONFID BAND IS
.05 SO TOTAL P= .1)

.1

t= 1.71713909197

REPLICATES ON UNKNOWN SAMPLE 3

y(c)= 0.0307393706507

x(d)= 0.00103307514615

ARMY CERTIFICATION
TETRYL
FOUND CONC




Figure I-9. Tetryl - Graph of Target-Found Concentration Points

Table I-21. DPA - Found Concentrations

ARMY CERTIFICATION	
DPA	FOUND CONC
(1)	5 4.31 (2) 5 4.3
(3)	10 7.5 (4) 10 7.6
(5)	20 15.93 (6) 20 16.96
(7)	50 35.01 (8) 50 35.1
(9)	100 79.74 (10) 100 81.62

Table I-22. DPA- Analysis of Target-Found Concentration Points

ARMY CERTIFICATION

DPA ANALYSIS OF 10 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 37 SD= 37.0585122925

FOUND CONC
MEAN= 28.807 SD= 29.5748620923

NO. RUNS 1 TOTAL X-Y ALL RUNS 10 NO. CONCENTR 10
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.633624595469
SLOPE= 0.795692556634

USE FOR ACCURACY

R= 0.997835330176

MEAN SQR DEV OF POINTS FROM REGRESSION= 5.82586021445

ST ERROR EST= 2.4136818793

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 8

TWO TAIL P LEVEL IS .1

t= 1.85954407943

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U U

REPLICATES= 1

Y(C)= 4.30510937847

X(D)= 12.2547388662

ARMY CERTIFICATION

DPA

FOUND CONC

81.52




Figure I-10. DPA -- Graph of Target-Found Concentration Points

Table I-23. 2,6-DNT - Found Concentrations -
HPLC Method

ARMY CERTIFICATION

26DNT

FOUND CONC

(1) 0 0 (2) 0 0

(3) 258 338 (4) 258 302

(5) 515 460 (6) 515 512

(7) 1030 1246 (8) 1030 1222

(9) 2575 2660 (10) 2575 2560

(11) 5150 5760 (12) 5150 5200

Table I-24. 2,6-DNT - Analysis of Target-Found
Concentration Points - HPLC Method

ARMY CERTIFICATION
26DNT
ANALYSIS OF 12 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 1588 SD= 1878.84635784

FOUND CONC
MEAN= 1695 SD= 1998.49698869

NO. RUNS 1 TOTAL X-Y ALL RUNS 12 NO. CONCENTR 12
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 9.88878232945
SLOPE= 1.86115315974
USE FOR ACCURACY

R= 0.997621596909
MEAN SQR DEV OF POINTS FROM REGRESSION= 20873.6483813
ST ERROR EST= 144.477155223
USE FOR PRECISION

T FOR CONFIDENCE BAND
D.F.= 10

TWO TAIL P LEVEL IS .1
t= 1.81245868646
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C
y(c)= 290.49849589
x(d)= 524.784244692 NG/2ML

ARMY CERTIFICATION
26DNT
FOUND CONC

5760.00




Figure I-11. 2,6-DNT - Graph of Target-Found Concentration Points - HPLC Method

Table I-25. 2,6-DNT-Found Concentrations -
GC/ECD Method

ARMY CERTIFICATION GC/ECD METHOD

2,6 DNT
FOUND CONC

(1)	0.84	0.03	(2)	0.04	0.06
(3)	0.84	0	(4)	0.04	0
(5)	0.87	0.03	(6)	0.07	0.03
(7)	0.07	0.88	(8)	0.07	0.04
(9)	0.14	0.16	(10)	0.14	0.21
(11)	0.14	0.12	(12)	0.14	0.06
(13)	0.36	0.45	(14)	0.36	0.34
(15)	0.36	0.39	(16)	0.36	0.41
(17)	0.72	0.84	(18)	0.72	0.73
(19)	0.72	0.71	(20)	0.72	0.71
(21)	0	0	(22)	0	0
(23)	0	0	(24)	0	0

Table I-26. 2,6-DNT - Analysis of Target-Found
Concentration Points

ARMY CERTIFICATION GC/ECD METHOD

2,6 DNT

ANALYSIS OF 24 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 0.22166666667 SD= 0.256831304246

FOUND CONC

MEAN= 0.225 SD= 0.277567069842

N0. RUNS 1 TOTAL X-Y ALL RUNS 24 N0. CONCENTR 24
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.0119009535528

SLOPE= 1.06872610625

R= 0.988886469517

MEAN SQR DEV OF POINTS FROM REGRESSION (ST ERROR EST)= 0.00178034051461
COMPUTE 1

D.F. = 22

ENTER 2 TAIL P LEVEL (USUALLY .1, EACH CONFID BAND IS
.05 SO TOTAL P= .1)

.1

t= 1.71713909197

REPLICATES ON UNKNOWN SAMPLE 3

y(c)= 0.0343436272298

x(d)= 0.0054574529401

ARMY CERTIFICATION GC/ECD METHOD
2,6-DNT
FOUND CONC




Figure I-12. 2,6-DNT - Graph of Target-Found
Concentration Points - GC/ECD Method

Table I-27. NG-Found Concentrations

ARMY CERTIFICATION	
NG	FOUND CONC
(1) 0 0 (2) 0 0	
(3) 7600 10120 (4) 7600 10740	
(5) 15360 13940 (6) 15360 17760	
(7) 30720 37600 (8) 30720 34260	
(9) 76000 78600 (10) 76000 78800	
(11) 153600 153980 (12) 153600 152560	

Table I-28. NG - Analysis of Target-Found
Concentration Points

ARMY CERTIFICATION

NG

ANALYSIS OF 12 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 47360 SD= 56038.9682599

FOUND CONC

MEAN= 49043.3333333 SD= 55472.3674591

NO. RUNS 1 TOTAL X-Y ALL RUNS 12 NO. CONCENTR 12
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 2199.55036995

SLOPE= 0.989100147031

USE FOR ACCURACY

R= 0.999202923622

MEAN SQR DEV OF POINTS FROM REGRESSION= 5393900.17082

ST ERROR EST= 2322.4771626

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 10

TWO TAIL P LEVEL IS .1

t= 1.81245868646

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

Y(C)= 6710.21417725

X(D)= 9875.83941875 NG/2ML

ARMY CERTIFICATION

NG
FOUND CONC
153980.00




Figure I-13. NG - Graph of Target-Found Concentration Points

Table I-29. PETN-Found Concentrations

ARMY CERTIFICATION	
PETN	FOUND CONC
(1)	1.2 1.5 (2) 1.2 1.06
(3)	1.2 1.73 (4) 1.2 2.32
(5)	2.4 2.59 (6) 2.4 2.07
(7)	2.4 3.52 (8) 2.4 3.42
(9)	4.8 5.24 (10) 4.8 5.42
(11)	4.8 5.98 (12) 4.8 7.07
(13)	12 12.22 (14) 12 11.6
(15)	12 13.2 (16) 12 12.56
(17)	24 25.03 (18) 24 24.19
(19)	24 26.41 (20) 24 27.34

Table I-30. PETH - Analysis of Target-Found
Concentration Points

ARMY CERTIFICATION
PETH
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 8.88 SD= 8.65682572188

FOUND CONC
MEAN= 9.7235 SD= 9.0930327724

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.437787216828
SLOPE= 1.04568837648
USE FOR ACCURACY
R= 0.99353506421
MEAN SQR DEV OF POINTS FROM REGRESSION= 0.779368060054
ST ERROR EST= 0.882818248596
USE FOR PRECISION
T FOR CONFIDENCE BAND
D.F.= 18
TWO TAIL P LEVEL IS .1
t= 1.73486096488
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U U
REPLICATES= 3
Y(C)= 1.36703361921
X(D)= 1.76926827162

ARMY CERTIFICATION

PETN
FOUND CONC

27.34




Figure I-14. PETN - Graph of Target-Found Concentration Points

TABLE II-1
QUANTITATIVE CERTIFICATION TESTING
STATISTICAL DATA SUMMARY

	Metal		Concrete		Brick		Transite	
	D.L. μg/cm ²	% Rec	D.L. μg/cm ²	% Rec	D.L. μg/cm ²	% Rec	D.L. μg/cm ²	% Rec
DNP	0.33	96%	1.74 (0.35)	31% 34%) ²	1.59	55%	2.22 (1.15)	34% 30%) ²
RDX	0.25	96%	0.63	78%	2.11	72%	3.48 (0.80)	67% 82%) ²
TNB	0.28	95%	0.98	75%	2.12	74%	3.46 (0.76)	53% 71%) ²
2,4-DNT	0.90	84%	1.08	78%	2.10	69%	3.52 (0.62)	65% 79%) ²
2,4,6-TNT	0.60	97%	1.57	74%	1.68	66%	3.18 (0.68)	51% 60%) ²
Tetryl	1.95	88%	4.29	51%	2.60	68%	4.13	56%
DPA	0.50	94%	2.44 (1.15)	74% 84%) ²	2.14	69%	3.71 (0.90)	67% 82%) ²
2,6-DNT	2.85	94%	6.46	86%	6.36	49%	2.04	79%
NG	9.36	94%	21.7	76%	32.5	44%	26.1	72%
PETN	10.2	86%	5.39	83%	20.6	62%	10.0	82%

¹ Calculated from four days of target vs. found concentrations using the procedures specified in the 1980 USATHAMA QA Plan.

² Calculated from three days of target vs. found concentrations.

QUANTITATIVE METHOD FOR THE DETERMINATION OF DNP,
RDX, TNB, 2,4-DNT, TNT ETRYL, AND DPA ON SURFACES

1. Application

Method used to extract the following compounds from metal, brick, concrete, transite surfaces:

2,4,-dinitrophenol DNP
cyclotrimethylenetrinitramine RDX
1,3,5-trinitrobenzene TNB
2,4-dinitrotoluene 2,4-DNT
2,4,6-trinitrotoluene 2,4,6-TNT
2,4,6-trinitrophenylmethylnitramine Tetryl
diphenylamine DPA

A. Tested Concentration Range

DNP	0.25 $\mu\text{g}/\text{cm}^2$ - 5.00 $\mu\text{g}/\text{cm}^2$
RDX	0.25 $\mu\text{g}/\text{cm}^2$ - 5.00 $\mu\text{g}/\text{cm}^2$
TNB	0.25 $\mu\text{g}/\text{cm}^2$ - 5.00 $\mu\text{g}/\text{cm}^2$
2,4-DNT	0.25 $\mu\text{g}/\text{cm}^2$ - 5.00 $\mu\text{g}/\text{cm}^2$
TNT	0.25 $\mu\text{g}/\text{cm}^2$ - 5.00 $\mu\text{g}/\text{cm}^2$
Tetryl	0.25 $\mu\text{g}/\text{cm}^2$ - 5.00 $\mu\text{g}/\text{cm}^2$
DPA	0.25 $\mu\text{g}/\text{cm}^2$ - 5.00 $\mu\text{g}/\text{cm}^2$

B. Sensitivity

Instrument response for each analyte is given below:

<u>Analyte</u>	<u>Concentration</u>	<u>Response</u>
DNP	25.04 ng/mL	3700 area units
RDX	24.86 ng/mL	1388 area units
TNB	24.97 ng/mL	4437 area units
2,4-DNT	24.99 ng/mL	6090 area units
TNT	24.91 ng/mL	4104 area units
Tetryl	25.30 ng/mL	2895 area units
DPA	24.91 ng/mL	2101 area units

C. Detection Limit

See Table II-1.

D. Interferences

Interferences present in some brick and transite samples were apparently random rather than systematic. For example, the HPLC analysis of extraction B of one blank (unspiked) brick surface (extracted Feb 9)

indicated the presence of a compound with a retention time of 1510 seconds. This compound interfered with the TNT. However, this interference was not observed in any other blank brick sample. Therefore, only the brick samples analyzed on Feb 9 were corrected for the interference.

E. Analysis Rate

Six samples can be extracted and prepared for analysis in three hours. Rate of analysis is given below, excluding calibration standards:

DNP	8 samples in an 8 hour day
RDX	8 samples in an 8 hour day
TNB	8 samples in an 8 hour day
2,4-DNT	8 samples in an 8 hour day
TNT	8 samples in an 8 hour day
Tetryl	8 samples in an 8 hour day
DPA	8 samples in an 8 hour day

2. Chemistry

2,4,-dinitrophenol C₆H₄N₂O₅
CAS RN 51-28-5
MP 112-114C

Cyclotrimethylenetrinitramine C₃H₆N₆O₆
CAS RN 121-82-4
MP 205-206C

1,3,5-trinitrobenzene C₆H₃N₃O₆
CAS RN 99-35-4
MP 122.5C

2,4-dinitrotoluene C₇H₆N₂O₄
CAS RN 121-14-2
MP 71C

2,4,6-trinitrotoluene C₇H₅N₃O₆
CAS RN 118-96-7
MP 80.1C

2,4,6-trinitrophenylmethylnitramine
CAS RN 479-45-8
MP 130C Explodes 187C

diphenylamine C₁₂H₁₁N
CAS RN 122-39-4
MP 53-54C BP 302C

3. Apparatus

A. Instrumentation

Waters Associates Model 6000A Solvent Delivery System
Waters Associates Model M-45 Solvent Delivery System
Waters Associates Model 660 Solvent Programmer
Waters Associates Model 440 Absorbance Detector
Waters Associates Intelligent Sample Processor (WISP)
Spectra-Physics Minigrator
Hewlett Packard 7133A Recorder

B. Parameters

Column: Spherisorb ODS 5 μ , 250 x 4.6 mm ID
Precolumn: Pellicular LC-18, 40 μ , 50 x 4.6 mm ID
Solvent System: linear gradient
 Initial: 30/70 CH₃CN/0.08 M acetic acid adjusted
 to pH 3.1 with NH₄OH
 Final: 50/50 CH₃CN/0.08 M acetic acid adjusted
 to pH 3.1 with NH₄OH
 Time: 35 minutes
Detector: UV at 254 nm
Flow Rate: 1.0 mL/min
Attenuation: 0.01 AUFS
Injection Volume: 100 μ L

C. Hardware/Glassware

Westinghouse Ultrasonic cleaner
8 ounce jars with teflon lined caps
25 mL graduated cylinders
microliter syringes
volumetric flasks - 50, 10, 5 mL
vials - WISP and 14 mL, with teflon lined caps

D. Chemicals

acetonitrile, HPLC grade
nitrogen
acetic acid
ammonium hydroxide
Standard Analytical Reference Material for each analyte

4. Standards

A. Calibration Standard

Stock Solution A: DNP, RDX, TNB, 2,4-DNT, TNT, TETRYL, DPA
Prepare individual stock solutions of 5.0 mg/mL. Combine
500 μ L individual stocks and dilute to 5 mL.

Stock Solution B:

Dilute 625 μ L Stock Solution A to 5 mL. Concentration is 62.5
 μ g/mL of each analyte

Calibration standards prepared in 25 mL volumetrics adding H₂O
so final solution is 55% H₂O/45% CH₃CN:

<u>Cal Std.</u>	<u>μL Stock B Added</u>	<u>Concentration Each Analyte</u>
1	10	25.0 μ g/L
2	20	50.0 μ g/L
3	40	100.0 μ g/L
4	80	200.0 μ g/L
5	200	500.0 μ g/L
6	400	1000.0 μ g/L

B. Control Spikes

Spiking stock solutions were prepared using 0.5 mg/mL stock solution.

<u>μL Stock Soln/x mL CH₃CN</u>	<u>Concentration Each Analyte</u>
125 μ L/10 mL	6.2 mg/L
250 μ L/10 mL	12.0 mg/L
500 μ L/10 mL	25.0 mg/L
500 μ L/5 mL	50.0 mg/L
625 μ L/5 mL	62.5 mg/L
1250 μ L/5 mL	125.0 mg/L
2500 μ L/5 mL	250.0 mg/L

400 μ L 6.2, 12.5, 25.0, 62.5, 125.0 mg/L Stocks spiked onto 10 cm²
concrete and metal surfaces.

200 μ L 12.5, 25.0, 50.0, 125.0, 250.0 mg/L stocks spiked onto 10 cm²
brick and transite surfaces.

Concentration of analytes on surface after spiking: 2.5 μ g, 5.0 μ g,
10.0 μ g, 25.0 μ g, 50 μ g.

5. Procedure

Extraction A

1. Spike 10 cm² surface sample with acetonitrile spike solution (volume dependent on surface type). Allow solvent to evaporate.
2. Transfer sample to 8 ounce jar and add 20 mL CH₃CN. Cover jar with teflon lined cap.
3. Sonicate for 10 minutes.
4. Transfer extract A to 50 mL volumetric flask, add 27 mL 0.08 M acetic acid and bring to volume with CH₃CN. Save surface sample for extraction B.

Extract A ready for analysis.

Extraction B

1. Add 20 mL CH₃CN to jar with surface sample. Sonicate for 10 minutes.
2. Transfer surface sample to a second jar. Add 20 mL CH₃CN and sonicate for 10 minutes more.
3. Combine extracts from steps 1 and 2 and evaporate using nitrogen to less than 5 mL.
4. Transfer evaporated extract to 10 mL volumetric flask. Add 5.5 mL to 0.08 M acetic acid and bring to volume with CH₃CN.

Extract B ready for analysis.

6. Calculations

Calculate found concentration for each analyte in each sample extract from daily calibration data.

Multiply found concentration by extract volume to find total µg in extract. Combine total µg in extracts A and B to find total µg on surface.

7. References

Lakings, D.B., Baker, R.J., and Crook, M.V., "Precision and Accuracy Assessment of the High Performance Liquid Chromatographic Analytical Technique for the Determination of Dinitrophenol (DNP); Cyclotrimethylene trinitramine (RDX); 1,3-Dinitrobenzene (DNB); 1,3,5-Trinitrobenzene (TNB); 2,4-Dinitrotoluene (2,4-DNT); Trinitrotoluene (TNT); 2,4,6-Trinitrophenyl-methylnitramine (Tetryl); and Diphenylamine (DPA)", Midwest Research Institute Technical Report No. 1, USATHAMA Contract No. DAAK11-81-C-0007, March, 1981.

QUANTITATIVE METHOD FOR THE DETERMINATION OF
2,6-DNT AND NG ON SURFACES

1. Application

Method used to extract the following compounds from metal, brick, concrete, transite surfaces:

2,6-dinitrotoluene 2,6-DNT
nitroglycerine NG

A. Tested Concentration Range

2,6-DNT	1.00 $\mu\text{g}/\text{cm}^2$ to 20.00 $\mu\text{g}/\text{cm}^2$
NG	12.50 $\mu\text{g}/\text{cm}^2$ to 125.00 $\mu\text{g}/\text{cm}^2$

B. Sensitivity

Instrument response for each analyte is given below:

<u>Analyte</u>	<u>Concentration</u>	<u>Response</u>
2,6-DNT	0.10 $\mu\text{g}/\text{mL}$	208750 area units
NG	1.25 $\mu\text{g}/\text{mL}$	186010 area units

C. Detection Limit

See Table II-1.

D. Interferences

No interferences were observed

E. Analysis Rate

Six samples can be extracted and prepared for analysis in three hours. Rate of analysis is given below excluding calibration standards:

2,6-DNT	16 Samples in an 8 hour day
NG	16 samples in an 8 hour day

2. Chemistry

2,6-dinitrotoluene C7H6N2O4
CAS RN 606-20-2
MP 66C

Nitroglycerine C3H5N3O9
CAS RN 55-63-0
MP Stable form 13.5C

3. Apparatus

A. Instrumentation

Beckman Model 110A Solvent Metering Pump
Waters Associates Model 450 Variable Wavelength Detector
Waters Associates Model U6K Injector
Hewlett Packard 3390A Integrator/Recorder

B. Parameters

Column: Spherisorb ODS, 5 μ , 250 x 4.6 mm ID
Precolumn: Pellicular LC-18, 40 μ , 50 x 4.6 mm ID
Solvent System: 35/65 CH₃CN/0.005 M t-butyl ammonium hydroxide,
pH 6.5. adjusted with 1N H₃PO₄
Detector: UV at 230 nm
Flow Rate: 1.0 mL/min
Attenuation: 0.01 AUFS
Injection Volume: 100 μ L

C. Hardware/Glassware

Westinghouse Ultrasonic cleaner
8 ounce jars with teflon lined caps
25 mL graduated cylinders
microliter syringes
volumetric flasks - 50, 10, 5 mL
vials - WISP and 14 mL, with teflon lined caps

D. Chemicals

Acetonitrile, HPLC grade
nitrogen
phosphoric acid
t-butyl ammonium hydroxide

Standard Analytical Reference Material for each analyte.

4. Standards

A. Calibration Standards

Prepare individual stock solutions:
5 mg/mL 2,6-DNT
50 mg/mL NG

Stock Solution A:

Combine 200 μ L 2,6-DNT stock and 250 μ L NG stock and dilute to 10 mL CH_3CN . Concentration is 0.1 mg/mL 2,6-DNT and 1.25 mg/mL NG.

Calibration standards prepared in 10 mL volumetric flasks adding 50% H_2O /50% CH_3CN .

<u>Cal Std.</u>	<u>μL Stock A added</u>	<u>2,6-DNT</u>	<u>NG</u>
1	10	0.1 μ g/mL	1.2 μ g/mL
2	20	0.2 μ g/mL	2.5 μ g/mL
3	40	0.4 μ g/mL	5.0 μ g/mL
4	80	0.8 μ g/mL	10.0 μ g/mL
5	200	2.0 μ g/mL	25.0 μ g/mL

B. Control Spikes

Spike Solutions prepared following chart below:

<u>Spike Solution</u>	<u>Amount Stock</u>	<u>Dilute with CH_3CN to</u>	<u>Concentration</u>
1	1 mL of 5 mg/mL 2,6-DNT	5 mL	1 mg/mL 2,6-DNT
2	0.5 mL of 5 mg/mL 2,6-DNT	5 mL	0.5 mg/mL 2,6-DNT
3	0.5 mL of 5 mg/mL 2,6-DNT 0.625 mL of 50 mg/mL NG	5 mL	0.5 mg/mL 2,6-DNT 6.25 mg/mL NG
4	0.250 mL of 5 mg/mL 2,6-DNT 0.312 mL of 50 mg/mL NG	5 mL	0.25 mg/mL 2,6-DNT 3.12 mg/mL NG
5	0.200 mL of 5 mg/mL 2,6-DNT 0.250 mL of 50 mg/mL NG	5 mL	0.20 mg/mL 2,6-DNT 2.50 mg/mL NG
6	0.200 mL of 5 mg/mL 2,6-DNT 0.250 mL of 50 mg/mL NG	10 mL	0.10 mg/mL 2,6-DNT 1.25 mg/mL NG
7	0.100 mL of 5 mg/mL 2,6-DNT 0.125 mL of 50 mg/mL NG	10 mL	0.05 mg/mL 2,6-DNT 0.62 mg/mL NG
8	0.025 mL of 5 mg/mL 2,6-DNT 0.031 mL of 50 mg/mL NG	5 mL	0.025 mg/mL 2,6-DNT 0.31 mg/mL NG

400 μ L of spike solutions 8, 7, 6, 4, 2 spiked onto 10 cm^2 concrete and metal surfaces

200 μ L of spike solutions 7, 6, 5, 3, 1 spiked onto 10 cm^2 brick and transite surfaces

Concentration on surface after spiking:

NG - 125 µg, 250 µg, 500 µg, 1250 µg

2,6-DNT - 10 µg, 20 µg, 40 µg, 100 µg, 200 µg.

5. Procedure

Extraction A

1. Spike 10 cm² surface sample with acetonitrile spike solution (volume dependent on surface type). Allow solvent to evaporate.
2. Transfer sample to 8 ounce jar and add 20 mL CH₃CN. Cover jar with teflon lined cap.
3. Sonicate for 10 minutes.
4. Transfer extract A to 50 mL volumetric flask, add 25 mL H₂O and bring to volume with CH₃CN.

Save surface sample for Extraction B.

Extraction B

1. Add 20 mL CH₃CN to jar with surface sample. Sonicate for 10 minutes.
2. Transfer surface sample to a second jar. Add 20 mL CH₃CN and sonicate for 10 minutes more.
3. Combine extracts from steps 1 and 2 and evaporate using nitrogen to less than 5 mL.
4. Transfer evaporated extract to 10 mL volumetric flask. Add 5.0 mL H₂O and bring to volume with CH₃CN. Extract B ready for analysis.

6. Calculations

Calculate found concentration for each analyte in each sample extract from daily calibration data.

Multiply found concentration by extract volume to find total µg in extract. Combine total µg in extracts A and B to find total µg on surface.

7. References

Lakings, D.B., Baker, R.J., and Crook, M.V., "Precision and Accuracy Assessment of the High Performance Liquid Chromatographic Analytical Technique for the Determination of Nitrobenzene (NB), 2,6-Dinitrotoluene (2,6-DNT), Nitroglycerin (NG), and Picric Acid (PA), Midwest Research Institute Technical Report No. 2, USATHAMA Contract No. DAAK11-81-C-0007, May, 1981.

QUANTITATIVE METHOD FOR THE DETERMINATION
OF PETN ON SURFACES

1. Application

Method used to extract pentaerythrite tetranitrate (PETN) from metal, brick, concrete, transite surfaces.

A. Tested Concentration Range:

PETN 5.0 $\mu\text{g}/\text{cm}^2$ to 100.0 $\mu\text{g}/\text{cm}^2$

B. Sensitivity

Instrument response for PETN is given below:

<u>Concentration</u>	<u>Response</u>
0.50 $\mu\text{g}/\text{mL}$	98025 area units

C. Detection Limit

See Table II-1.

D. Interferences

There were no interferences.

E. Analysis Rate

Six samples can be extracted and prepared for analysis in three hours. Rate of analysis is given below, excluding calibration standards:

PETN 32 samples in an 8 hour day

2. Chemistry

Pentaerythrite tetranitrate C5H8N4O12
CAS RN 78-11-5
MP 140-141 C

3. Apparatus

A. Instrumentation

Beckman Model 110A Solvent Metering Pump
Waters Associates Model 450 Variable Wavelength Detector
Waters Associates Model U6K Injector
Hewlett Packard 3390A Integrator/Recorder

B. Parameters

Column: Spherisorb ODS, 5 μ , 250 x 4.6 mm ID
Precolumn: Pellicular LC-18, 40 μ , 50 x 4.6 mm ID
Solvent System: 65% CH₃CN/35% H₂O
Detector: UV at 230 nm
Flow Rate: 1.0 mL/min
Attenuation: 0.01 AUFS
Injection Volume: 100 μ L

C. Hardware/Glassware

Westinghouse Ultrasonic cleaner
8 ounce jars with teflon lined caps
25 mL graduated cylinders
microliter syringes
volumetric flasks - 50, 10, 5 mL
vials - WISP and 14 mL, with teflon lined caps

D. Chemicals

Acetonitrile, HPLC grade
Standard Analytical Reference Material for PETN

4. Standards

A. Calibration Standards:

Prepare stock solution as follows:
200 μ L of SARM (50 mg/mL) in 10 mL CH₃CN = 1.0 mg/mL

Calibration Standards prepared in 10 mL volumetric flasks adding
50% H₂O/50% CH₃CN.

<u>Cal. Std</u>	<u>μL Stock added</u>	<u>Concentration PETN</u>
1	5	0.5 μ g/mL
2	10	1.0 μ g/mL
3	20	2.0 μ g/mL
4	40	4.0 μ g/mL
5	100	10.0 μ g/mL
6	200	20.0 μ g/mL

B. Control Spikes

Spike solutions prepared following chart below:

<u>Spike Solution</u>	<u>Amount Stock</u>	<u>Dilute with CH₃CN to</u>	<u>Concentration</u>
1	12.5 μ L of 50 mg/mL	5 mL	0.125 mg/mL
2	50 μ L of 50 mg/mL	10 mL	0.25 mg/mL
3	100 μ L of 50 mg/mL	10 mL	0.50 mg/mL
4	100 μ L of 50 mg/mL	5 mL	1.0 mg/mL
5	125 μ L of 50 mg/mL	5 mL	1.25 mg/mL
6	500 μ L of 50 mg/mL	10 mL	2.5 mg/mL
7	500 μ L of 50 mg/mL	5 mL	5.0 mg/mL

400 μ L spike solutions 1,2,3,5 and 6 spiked onto 10 cm² metal surfaces

200 μ L spike solutions 2,3,4,6,7 spiked onto 10 cm² concrete, transite, and brick surfaces.

Concentration of analytes on surface after spiking: 50 μ g, 100 μ g, 200 μ g, 500 μ g, 1000 μ g.

5. Procedure

1. Spike 10 cm² surface sample with acetonitrile spike solution (volume dependent on surface type). Allow solvent to evaporate.
2. Transfer sample to 8 ounce jar.
3. Add 12 mL CH₃CN. Cover jar with teflon lined cap. Sonicate for 10 min.
4. Transfer extract to 50 mL volumetric flask.
5. Repeat steps 3 and 4 twice, adding extracts to same 50 mL volumetric flask.
6. Rinse jar with 17 mL H₂O and transfer to volumetric flask.
7. Bring to volume with CH₃CN.
8. Ready for HPLC analysis.

6. Calculations

Calculate found concentration for each analyte in each sample extract from daily calibration data.

Multiply found concentration by extract volume to find total μ g in extract.

7. References

None

Table II-2. DNP on Metal - Target US. Found Concentrations

2,4-DINITROPHENOL (24DNP)	
METAL SURFACE	
TARGET CONC. us FOUND CONC	us FOUND CONC
us/10 sq cm	us/10 sq cm
2.500	1.920
	2.270
	2.260
	2.470
5.000	3.700
	4.030
	4.610
	5.100
10.000	9.060
	8.020
	10.040
	11.220
25.000	23.500
	22.460
	25.240
	24.610
50.000	46.870
	47.850
	48.840
	48.220

Table II-3. DNP on Metal - Analysis of Target-
Found Concentration Points

2,4-DINITROPHENOL (24DNP)
METAL SURFACE
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.0350535872

FOUND CONC
MEAN= 17.6145 SD= 17.4133452982

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.226280339806
SLOPE= 0.964366504854
USE FOR ACCURACY
R= 0.998797261237
MEAN SQR DEV OF POINTS FROM REGRESSION= 0.769459159271
ST ERROR EST= 0.877188212
USE FOR PRECISION
T FOR CONFIDENCE BAND
D.F.= 18
TWO TAIL P LEVEL IS .1
t= 1.73406096408
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C
(EACH TARGET CONC CONSIDERED INDEP SAMPLE
MEASURED 1 TIME(S))
Y(C)= 1.37295729236
X(D)= 3.30309425069

2,4-DINITROPHENOL (24DNP)
METAL SURFACE
FOUND CONC




Figure 11-1. DNP on Metal - Graph of Target-Found Concentration Points

Table II-4. DYP on Metal - Inaccuracy and Imprecision Data

**2,4-DINITROPHENOL (24DNP)
METAL SURFACE
STATISTICAL DATA USED TO DETERMINE PERCENT
INACCURACY AND IMPRECISION**

Target Conc ug/10 sq cm	Mean Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	2.230	0.220	-10.800	10.232
5.000	4.360	0.620	-12.900	14.229
10.000	9.585	1.367	-4.150	14.261
25.000	23.953	1.220	-4.190	5.126
50.000	47.945	0.825	-4.110	1.720
Means		0.854	-7.210	9.114

2,4-DINITRUPHENOL (24DNP)

METAL SURFACE

MEAN INACCURACY

CONCENTRATION (REPORT)

-4.11

-12.80

TARGET CONC

INTERVAL = 0.869

Graph of Inaccuracy

Figure II-2. DNP on Metal

Arthur D. Little, Inc.

2,4-DINITROPHENOL (24DNP)
 METAL SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-3. DNP on Metal - Graph of Imprecision

Table II-5. DNP on Concrete - Target U.S. Found Concentrations

2,4-DINITROPHENOL (24DNP) CONCRETE SURFACE TARGET CONC. US FOUND CONC	
Target Conc ug/10 sq cm	Found Conc ug/10 sq cm
2.500	0.710 1.010 1.470 1.250
5.000	1.320 1.890 2.690 1.920
16.000	2.270 3.500 3.800 3.290
25.000	6.210 8.430 8.980 8.560
50.000	10.800 17.570 17.540 16.980

Table II-6. DNP on Concrete - Analysis of Target-
Found Concentration Points

2,4-DINITROPHENOL (24DNP)
CONCRETE SURFACE
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.0350535072

FOUND CONC
MEAN= 6.0095 SD= 5.74332424029

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.308970469256
SLOPE= 0.308136731392
USE FOR ACCURACY
R= 0.967603816639
MEAN SQR DEV OF POINTS FROM REGRESSION= 2.219410059
ST ERROR EST= 1.48977141166
USE FOR PRECISION
T FOR CONFIDENCE BAND
D.F.= 18
TWO TAIL P LEVEL IS .1
t= 1.73406096408
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C
Y(C)= 3.02503353561
X(D)= 17.4060970424

2,4-DINITROPHENOL (24DNP)
CONCRETE SURFACE
FOUND CONC




Figure II-4. DNP on Concrete - Graph of Target-Found Concentration Points

Table II-7. DNP on Concrete - Inaccuracy and Precision Data

2,4-DINITROPHENOL (24DNP) CONCRETE SURFACE STATISTICAL DATA USED TO DETERMINE PERCENT INACCURACY AND IMPRECISION					
Mn Target Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision	
2.500	1.110	0.326	-55.600	29.386	
5.000	1.955	0.562	-60.900	20.768	
10.000	3.215	0.664	-67.830	20.649	
25.000	8.045	1.246	-67.820	15.483	
50.000	15.723	3.293	-68.555	20.944	
Means		1.218	-64.145	23.046	

2,4-DINITROPHENOL (24DNP)
 CONCRETE SURFACE
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure II-5. DNP on Concrete - Graph of Inaccuracy

2,4-DINITROPHENOL (24DNP)
 CONCRETE SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-6. DNP on Concrete - Graph of Imprecision

Table II-8. DNP on Concrete (3 days) - Target vs.
Found Concentrations

2,4-DINITROPHENOL (24DNP) CONCRETE SURFACE	
TARGET CONC. US FOUND CONC	
Target Conc ug/10 sq cm	Found Conc ug/10 sq cm
2.500	1.010 1.470 1.250
5.000	1.890 2.690 1.920
10.000	3.500 3.800 3.290
25.000	8.430 8.980 8.560
50.000	17.570 17.540 16.980

Table II-9. DNP on Concrete (3 days) - Analysis of Target-Found Concentration Points

2,4-DINITROPHENOL (24DNP)
CONCRETE SURFACE

ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.1953683274

FOUND CONC
MEAN= 6.592 SD= 6.17940727625

NØ. RUN: 1 TOTAL X-Y ALL RUNS 15 NØ. CONCENTR 15
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.31678802589

SLOPE= 0.339200647249

USE FOR ACCURACY

R= 0.998781992786

MEAN SQR DEV OF POINTS FROM REGRESSION= 0.100113722929

ST ERROR EST= 0.316407326662

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F. = 13

TWO TAIL P LEVEL IS .1

t= 1.77093170942

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 0.915196314813

x(d)= 3.50862633373

2,4-DINITROPHENOL (24DNP)
CONCRETE SURFACE
FOUND CONC

17.57

$y(c)$

$y(d)$

3

0.00

VERTICAL AXIS TIC INTERVAL= 1.757

TARGET CONC

50

40

30

20

10

0

Figure II-7. DNP on Concrete (3 days) - Graph of Target-Found Concentration Points

Table II-10. DNP on Concrete (3 days) - Inaccuracy and Imprecision Data

2,4-DINITROPHENOL (24DNP) CONCRETE SURFACE STATISTICAL DATA USED TO DETERMINE PERCENT INACCURACY AND IMPRECISION					
Mn Target Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision	
2.500	1.243	0.230	-50.267	18.504	
5.000	2.167	0.453	-56.667	20.929	
10.000	3.530	0.256	-64.701	7.261	
25.000	8.657	0.287	-65.373	3.321	
50.000	17.363	0.332	-65.273	1.914	
Means		0.312	-60.456	10.386	

2,4-DINITROPHENOL (24DNP)
 CONCRETE SURFACE
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure II-8. DNP on Concrete (3 days) - Graph of Inaccuracy

2,4-DINITROPHENOL (24DNP)
 CONCRETE SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-9. DNP on Concrete (3 days) - Graph of Imprecision

Table II-11. DNP on Brick - Target vs. Found Concentrations

2,4-DINITROPHENOL (24DNP) BRICK SURFACE	
TARGET CONC. ug/10 sq cm	VS. FOUND CONC. Found Conc ug/10 sq cm
2.500	1.000 1.360 0.810 1.110
5.000	1.690 4.000 2.780 3.620
10.000	5.560 3.500 4.890 6.800
25.000	13.450 10.200 11.320 14.100
50.000	28.950 34.090 21.730 26.490

Table II-12. DNP on Brick - Analysis of Target-Found Concentration Points

2,4-DINITROPHENOL (24DNP)
BRICK SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.0350535872

FOUND CONC
MEAN= 9.8725 SD= 18.2455330148

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.352117718446
SLOPE= 0.532682038835

USE FOR ACCURACY
R= 0.972873761094

MEAN SQR DEV OF POINTS FROM REGRESSION= 5.92937890036
ST ERROR EST= 2.43503160151

USE FOR PRECISION
T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

Y(C)= 4.08728778129

X(d)= 15.86533395681

Table II-13. DNP on Brick - Inaccuracy and Imprecision Data

2,4-DINITROPHENOL (24DNP)

BRICK SURFACE
STATISTICAL DATA USED TO DETERMINE PERCENT
INACCURACY AND IMPRECISION

Mn Target Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	2.500	0.230	-57.200	21.461
5.000	3.023	1.024	-39.550	33.885
10.000	5.188	1.375	-48.125	26.514
25.000	12.268	1.819	-50.930	14.830
50.000	27.815	5.146	-44.370	18.501
Means		1.919	-48.035	23.030

2,4-DINITROPHENOL (24DNP)
 BRICK SURFACE
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure II-11. DNP on Brick - Graph of Inaccuracy

2,4-DINITROPHENOL (24DNP)
 BRICK SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-12. DNP on Brick - Graph of Imprecision

Table II-14. DNP on Transite - Target vs. Found Concentrations

2,4-DINITROPHENOL (24DNP)	
TRANSITE SURFACE	
TARGET CONC. VS. FOUND CONC	
Target Conc ug/10 sq cm	Found Conc ug/10 sq cm
2.500	0.980
	0.710
	0.250
	1.280
5.000	2.260
	1.440
	1.200
	1.970
10.000	3.880
	2.380
	2.790
	2.630
25.000	9.700
	7.140
	6.960
	7.050
50.000	24.050
	16.740
	12.550
	15.680

Table II-15. DNP on Transite - Analysis of Target-
Found Concentration Points

2,4-DINITROPHENOL (24DNP)
TRANSITE SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.0350535872

FOUND CONC
MEAN= 6.002 SD= 6.55535747148

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.297536407767

SLOPE= 0.344839805825

USE FOR ACCURACY

R= 0.948720859866

MEAN SQR DEV OF POINTS FROM REGRESSION= 4.53277563377

ST ERROR EST= 2.12903161078

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(d) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

Y(c)= 3.58398803391

X(d)= 22.2382830955

Table II-16. DNP on Transite - Inaccuracy and Imprecision Data

2,4-DINITROPHENOL (24DNP)

TRANSITE SURFACE

STATISTICAL DATA USED TO DETERMINE PERCENT

INACCURACY AND IMPRECISION

Mn Target Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	0.805	0.437	-67.800	54.304
5.000	1.718	0.484	-65.650	28.183
10.000	2.920	0.662	-70.800	22.667
25.000	7.713	1.327	-69.150	17.206
50.000	17.235	4.867	-65.490	28.205
Means		1.555	-67.770	30.113

2,4-DINITROPHENOL (24DNP)
 TRANSITE SURFACE
 MEAN INACCURACY

MEAN FOUND CONCENTRATION (REPORT)
 -65.49




Figure II-14. DNP on Transite - Graph of Inaccuracy

2,4-DINITROPHENOL (24DNP)
 TRANSITE SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-15. DNP on Transite - Graph of Imprecision

Table II-17. DNP on Transite (3 days) - Target vs. Found Concentrations

2,4-DINITROPHENOL (24DNP)	
TRANSITE SURFACE	
TARGET CONC. US FOUND CONC	
ug/10 sq cm	ug/10 sq cm
2.500	0.710 0.250 1.280
5.000	1.440 1.200 1.970
10.000	2.380 2.790 2.630
25.000	7.140 6.960 7.030
50.000	16.740 12.550 15.680

Table II-18. DNP on Transite (3 days) - Analysis of Target-Found Concentration Points

2,4-DINITROPHENOL (24DNP)
TRANSITE SURFACE

ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.1953683274

FOUND CONC
MEAN= 5.3846666667 SD= 5.52760328671

NO. RUNS 1 TOTAL X-Y ALL RUNS 15 NO. CONCENTR 15
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.161960623674

SLOPE= 0.299817691478

USE FOR ACCURACY

R= 0.986918388416

MEAN SQ DEV OF POINTS FROM REGRESSION= 0.855263021743

ST ERROR EST= 0.924804315378

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 13

TWO TAIL P LEVEL IS .1

t= 1.77093170942

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

Y(C)= 1.5870829987

X(D)= 11.5032476202

Table II-19. DNP on Transite (3 days) - Inaccuracy and Imprecision Data

2,4-DINITROPHENOL (24DNP) TRANSITE SURFACE STATISTICAL DATA USED TO DETERMINE PERCENT INACCURACY AND IMPRECISION				
Mn Target Conc ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	0.747	0.516	-70.133	69.104
5.000	1.537	0.394	-69.267	25.640
10.000	2.600	0.207	-74.000	7.948
25.000	7.050	0.090	-71.800	1.277
50.000	14.990	2.179	-70.020	14.533
Means		0.677	-71.044	23.700

ACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure II-17. DNP on Transit (3 days) - Graph of Inaccuracy

2,4-DINITROPHENOL (24DNP)
 TRANSITE SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-18. DNP on Transite (3 days) - Graph of Imprecision

Table II-20. RDX on Metal - Target vs. Found Concentrations

CYCLOTRIMETHYLENETRINITRAMINE (RDX)

METAL SURFACE	
TARGET CONC. ug/10 sq cm	VS. FOUND CONC. Found Conc ug/10 sq cm
2.500	2.340
	2.040
	2.180
	2.620
5.000	4.100
	4.250
	4.840
	5.070
25.000	24.180
	23.960
	23.070
	24.460
50.000	49.060
	46.620
	46.960
	48.510

Table II-21. RDX on Metal - Analysis of Target-Found Concentrations Points

CYCLOTRIMETHYLENETRINITRAMINE (RDX)
METAL SURFACE

ANALYSIS OF 16 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 20.625 SD= 19.6956001855

FOUND CONC
MEAN= 19.64125 SD= 18.9047401022

N0. RUNS 1 TOTAL X-Y ALL RUNS 16 N0. CONCENTR 16
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.144172932331

SLOPE= 0.959293233083

USE FOR ACCURACY

R= 0.999424264881

MEAN SQR DEV OF POINTS FROM REGRESSION= 0.440790601532
ST ERROR EST= 0.66392062894

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F. = 14

TWO TAIL P LEVEL IS .1

t= 1.7613101065

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C
(EACH TARGET CONC CONSIDERED INDEP SAMPLE

MEASURED 1 TIME(S))

Y(C)= 1.10196465687

X(D)= 2.58815970399

CYCLOTRIMETHYLENETRINITRAMINE (RDX)
METAL SURFACE
FOUND CONC




Figure II-19. RDX on Metal - Graph of Target-Found Concentration Points

Table II-22. RDX on Metal - Inaccuracy and Imprecision Data

CYCLOTRIMETHYLENETRINITRAMINE (RDX)

METAL SURFACE

STATISTICAL DATA USED TO DETERMINE PERCENT

INACCURACY AND IMPRECISION

Mn Target Conc ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	2.295	0.249	-8.200	10.847
5.000	4.565	0.464	-8.700	10.166
25.000	23.918	0.601	-4.330	2.512
50.000	47.788	1.182	-4.425	2.473
Means		0.624	-6.414	6.499




Figure 11-20. RDX on Metal - Graph of Inaccuracy

CYCLOTRIMETHYLENETRINITRAMINE (RDX)
 METAL SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-21. RDX on Metal - Graph of Imprecision

Table II-23. RDX on Concrete - Target vs. Found Concentrations

CYCLOTRIMETHYLENE TRINITRAMINE (RDX)

CONCRETE SURFACE	
TARGET CONC. vs. FOUND CONC.	
Target Conc ug/10 sq cm	Found Conc ug/10 sq cm
2.500	1.610
	1.740
	1.950
	2.090
5.000	2.860
	3.710
	4.320
	4.210
10.000	6.730
	7.650
	7.520
	7.930
25.000	19.370
	20.410
	20.410
	17.490
50.000	37.810
	41.120
	41.240
	35.160

Table II-24. RDX on Concrete - Analysis of Target-
Found Concentration Points

CYCLOTRIMETHYLENETRINITRAMINE (RDX)
CONCRETE SURFACE
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.0350535872

FOUND CONC
MEAN= 14.2665 SD= 14.1353562966

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.168798543689
SLOPE= 0.780286407767
USE FOR ACCURACY
R= 0.995533764773
MEAN SQR DEV OF POINTS FROM REGRESSION= 1.8713384477
ST ERROR EST= 1.36796580648
USE FOR PRECISION
T FOR CONFIDENCE BAND
D.F.= 18
TWO TAIL P LEVEL IS .1
t= 1.73406096408
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C
y(c)= 2.32519575231
x(d)= 6.34666977015

CYCLOTRIMETHYLENETRINITRAMINE (RDX)
 CONCRETE SURFACE
 FOUND CONC




Figure II-22. RDX on Concrete - Graph of Target-Found Concentration Points

Table II-25. Inaccuracy and Imprecision Data

**CYCLOTRIMETHYLENETRINITRAMINE (RDX)
CONCRETE SURFACE
STATISTICAL DATA USED TO DETERMINE PERCENT
INACCURACY AND IMPRECISION**

Mn Target Conc ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	1.848	0.214	-26.100	11.570
5.000	3.775	0.665	-24.500	17.623
10.000	7.450	0.514	-25.425	6.896
25.000	19.420	1.377	-22.320	7.090
50.000	38.833	2.919	-22.335	7.517
Means		1.138	-24.136	10.141




Figure 11-23. RDX on Concrete - Graph of Inaccuracy

CYCLOTRIMETHYLENENITRAMINE (RDX)
 CONCRETE SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-24. RDX on Concrete - Graph of Imprecision

Table II-26. RDX on Brick - Target vs. Found Concentrations

CYCLOTRIMETHYLENETRINITRAMINE (RDX)

CYCLOTRIMETHYLENETRINITRAMINE BRICK SURFACE TARGET CONC. VS. FOUND CONC.	
Target Conc ug/10 sq cm	Found Conc ug/10 sq cm
2.500	1.340 2.070 1.020 1.160
5.000	2.450 7.380 4.120 5.010
10.000	9.080 6.410 7.190 9.900
25.000	21.180 14.890 14.480 20.990
50.000	36.030 45.030 24.020 41.020

Table II-27. RDX on Brick - Analysis of Target-Found
Concentration Points

CYCLOTRIMETHYLENETRINITRAMINE (RDX)
BRICK SURFACE
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.0350535872

FOUND CONC
MEAN= 13.7385 SD= 13.6021253274

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.432779126214

SLOPE= 0.71922815534

USE FOR ACCURACY

R= 0.95362437934

MEAN SQR DEV OF POINTS FROM REGRESSION= 17.6939762945

ST ERROR EST= 4.20642084135

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(0) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U CB

y(c)= 8.10167664923

x(d)= 21.0600746217

CYCLOTRIMETHYLENETRINITRAMINE (RDX)
BRICK SURFACE
FOUND CONC




Figure 14-25. RDX on Brick - Graph of Target-Found Concentration Points

Table 11-28. RDX on Brick - Inaccuracy and Imprecision Data

CYCLOTRIMETHYLENETRINITRAMINE (RDX)

BRICK SURFACE
STATISTICAL DATA USED TO DETERMINE PERCENT
INACCURACY AND IMPRECISION

Mn Target Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	1.398	0.467	-44.100	33.422
5.000	4.740	2.055	-5.200	43.358
10.000	8.145	1.620	-18.550	19.894
25.000	17.885	3.700	-28.460	20.686
50.000	36.525	9.113	-26.950	24.951
Means		3.391	-24.652	28.462

CYCLOTRIMETHYLENETRINITRAMINE (RDX)
 BRICK SURFACE
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure 11-26. RDX on Brick - Graph of Inaccuracy

CYCLOTRIMETHYLENETRINITRAMINE (RDX)
 BRICK SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-27. RDX on Brick - Graph of Imprecision

Table II-29. RDX on Transite - Target vs. Found Concentrations

CYCLOTRIMETHYLENETRINITRAMINE (RDX)

TRANSITE SURFACE	
TARGET CONC. US FOUND CONC	
Target Conc Found Conc	
ug/10 sq cm ug/10 sq cm	
2.500	0.710
	1.640
	1.570
	1.410
5.000	1.610
	3.900
	3.670
10.000	3.480
	2.630
	7.590
25.000	7.150
	7.840
	6.640
	15.830
50.000	19.870
	18.270
	12.540
	39.670
	44.350
	38.070

Table 11-30. RDX on Transite - Analysis of Target-
Found Concentration Points

CYCLOTRIMETHYLENETRINITRAMINE (RDX)
TRANSITE SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.3 SD= 18.0350535872

FOUND CONC
MEAN= 11.922 SD= 13.6804011396

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.544051375404

SLOPE= 0.673840614887

USE FOR ACCURACY

R= 0.888327791599

MEAN SQR DEV OF POINTS FROM REGRESSION= 41.6587590559
ST ERROR EST= 6.45435969372

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 11.2231552201

x(d)= 34.829209243




Figure 11-28. RDX on Transit - Graph of Target-Found Concentration Points


Table II-31. RDX on Transite - Inaccuracy and Imprecision Data

CYCLOTRIMETHYLENETRINITRAMINE (RDX)
 TRANSITE SURFACE
 STATISTICAL DATA USED TO DETERMINE PERCENT
 INACCURACY AND IMPRECISION

Mn Target Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	1.333	0.426	-46.700	31.971
5.000	3.165	1.051	-36.700	33.200
10.000	6.303	2.465	-36.975	39.110
25.000	15.153	5.913	-39.390	39.024
50.000	33.658	14.328	-32.685	42.571
Means		4.837	-38.490	37.175

CYCLOTRIMETHYLENETRINITRAMINE (RDX)
 TRANSITE SURFACE
 MEAN FOUND CONCENTRATION (REPORT)
 MEAN INACCURACY

-32.69



-46.70

TARGET CONC

VERTICAL AXIS TIC INTERVAL = 1.4015

Figure 11-29. RDX on Transite - Graph of Inaccuracy

CYCLOTRIMETHYLENETRINITRAMINE (RDX)
 TRANSITE SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-30. RDX on Transite - Graph of Imprecision

Table II-32. RDX on Transite (3 days) - Target vs.
Found Concentrations

CYCLOTRIMETHYLENETRINITRAMINE (RDX)

TRANSITE SURFACE	
TARGET CONC. VS. FOUND CONC.	
Target Conc ug/10 sq cm	Found Conc ug/10 sq cm
2.500	1.640 1.570 1.410
5.000	3.900 3.670 3.480
10.000	7.590 7.150 7.840
25.000	15.830 19.870 18.270
50.000	39.670 44.350 38.070

Table II-33. RDX on Transite (3 days) - Analysis of Target-Found Concentration Points

CYCLOTRIMETHYLENETRINITRAMINE (RDX)
TRANSITE SURFACE

ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.1953683274

FOUND CONC
MEAN= 14.287333333 SD= 14.9412293751

N0. RUNS 1 TOTAL X-Y ALL RUNS 15 N0. CONCENTR 15
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.806491370011
SLOPE= 0.815882416397

USE FOR ACCURACY

R= 0.993578286336

MEAN SQR DEV OF POINTS FROM REGRESSION= 3.07780840387
ST ERROR EST= 1.75436837747

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F. = 13

TWO TAIL P LEVEL IS .1

t= 1.77093170942

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

Y(C)= 2.51147198116

X(D)= 8.04284789281

CYCLOTRIMETHYLENETRINITRAMINE (RDX)
TRANSITE SURFACE
FOUND CONC




Figure 11-31. RDX on Transite (3 days) - Graph of Target-Found Concentration Points

Table II-34. RDX on Transite (3 days) - Inaccuracy and Imprecision Data

CYCLOTRIMETHYLENETRINITRAMINE (RDX)

TRANSITE SURFACE
STATISTICAL DATA USED TO DETERMINE PERCENT
INACCURACY AND IMPRECISION

Mn Target Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	1.540	0.118	-38.400	7.656
5.000	3.683	0.210	-26.333	5.710
10.000	7.527	0.349	-24.733	4.641
25.000	17.990	2.035	-28.040	11.309
50.000	40.697	3.263	-18.607	8.019
Means		1.195	-27.223	7.467

CYCLOTRIMETHYLENETRINITRAMINE (RDX)
 TRANSITE SURFACE
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure 11-32. RDX on Transite (3 days) - Graph of Inaccuracy

CYCLOTRIMETHYLENETRINITRAMINE (RDX)
 TRANSITE SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-33. RDX on Transite (3 days) - Graph of Imprecision

Table II-35. TNB on Metal - Target vs. Found Concentrations

1,3,5-TRINITROBENZENE (135TNB)

METAL SURFACE	
TARGET CONC. ug/10 sq cm	VS FOUND CONC Found Conc ug/10 sq cm
2.500	2.130
	2.290
	2.280
	2.470
5.000	3.690
	4.460
	4.740
	4.820
10.000	9.380
	8.150
	9.600
	10.190
25.000	23.560
	23.480
	25.550
	24.340
50.000	46.760
	47.630
	46.340
	48.390

Table II-36. TNB on Metal - Analysis of Target-Found
Concentration Points

1,3,5-TRINITROBENZENE (135TNB)
METAL SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.0350535872

FOUND CONC
MEAN= 17.5125 SD= 17.1872997405

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.102313915858
SLOPE= 0.95215210356

USE FOR ACCURACY
R= 0.999116467986

MEAN SQR DEV OF POINTS FROM REGRESSION= 0.550752890156
ST ERROR EST= 0.742127273556

USE FOR PRECISION
T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 1.25068860685

x(d)= 2.83188868416

1,3,5-TRINITROBENZENE (135TNB)
METAL SURFACE
FOUND CONC

48.39




Figure II-34. TNB on Metal - Graph of Target-Found Concentration Points

Table II-37. TNB on Metal - Inaccuracy and Imprecision Data

1,3,5-TRINITROBENZENE (135TNB)

METAL SURFACE

STATISTICAL DATA USED TO DETERMINE PERCENT INACCURACY AND IMPRECISION

Mn Target Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	2.293	0.139	-8.300	6.069
5.000	4.428	0.515	-11.450	11.639
10.000	9.330	0.858	-6.700	9.194
25.000	24.233	0.960	-3.070	3.962
50.000	47.280	0.914	-5.440	1.934
Means		0.677	-6.992	6.560

1,3,5-TRINITROBENZENE (135TNB)
METAL SURFACE

MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)

-3.07




Figure 11-35. TNB on Metal - Graph of Inaccuracy

1,3,5-TRINITROBENZENE (135TNB)
 METAL SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION <REPORT>




Figure II-36. TNB on Metal - Graph of Imprecision

Table II-38. TNB on Concrete - Target vs. Found Concentrations

1,3,5-TRINITROBENZENE (135TNB)
CONCRETE SURFACE

TARGET CONC. US FOUND CONC	Target Conc	Found Conc
ug/10 sq cm	ug/10 sq cm	ug/10 sq cm
2.500	1.950	1.920
		2.090
		1.800
5.000	3.160	3.510
		5.120
		3.690
10.000	5.880	7.170
		8.570
		7.620
25.000	20.890	22.660
		20.540
		17.150
50.000	38.570	42.370
		34.610
		34.000

Table II-39. TNB on Concrete - Analysis of Target-
Found Concentration Points

1,3,5-TRINITROBENZENE (135TNB)
CONCRETE SURFACE
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.0350535072

FOUND CONC
MEAN= 14.1635 SD= 13.7366717317

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.222028721683
SLOPE= 0.753593042071
USE FOR ACCURACY
R= 0.989402029999
MEAN SQR DEV OF POINTS FROM REGRESSION= 4.19942063895
ST ERROR EST= 2.04924879869
USE FOR PRECISION
T FOR CONFIDENCE BAND
D.F.= 18

THO TAIL P LEVEL IS .1
t= 1.73406096408
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C
y(c)= 3.95809784769
x(d)= 9.81755206131

1,3,5-TRINITROBENZENE (135TNB)
CONCRETE SURFACE
FOUND CONC




Figure 11-37. TNB on Concrete - Graph of Target-Found Concentration Points

Table II-40. TNB on Concrete - Inaccuracy and Imprecision Data

1,3,5-TRINITROBENZENE (135TNB)					
CONCRETE SURFACE					
STATISTICAL DATA USED TO DETERMINE PERCENT					
INACCURACY AND IMPRECISION					
Mn Target Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision	
2.500	1.940	0.119	-22.400	6.142	
5.000	3.870	0.062	-22.600	22.271	
10.000	7.310	1.118	-26.900	15.291	
25.000	20.310	2.302	-18.760	11.334	
50.000	37.388	3.891	-25.225	10.406	
Means		1.658	-23.177	13.089	

1,3,5-TRINITROBENZENE (135TNB)
 CONCRETE SURFACE
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure 11-38. TNB on Concrete - Graph of Inaccuracy

1,3,5-TRINITROBENZENE (135TNB)
 CONCRETE SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-39. TNB on Concrete - Graph of Imprecision

Table II-41. TNB on Brick - Target vs. Found Concentrations

1,3,5-TRINITROBENZENE (135TNB)

BRICK SURFACE

TARGET CONC. VS. FOUND CONC.

Target Conc ug/10 sq cm	Found Conc ug/10 sq cm
2.500	1.390 2.180 1.270 2.340
5.000	2.800 6.180 4.090 4.950
10.000	8.250 6.340 7.600 9.480
25.000	22.460 14.690 12.960 21.150
50.000	39.220 43.160 23.840 43.340

Table II-42. TNB on Brick - Analysis of Target-Found
(Concentration Points)

1,3,5-TRINITROBENZENE (135TNB)

BRICK SURFACE
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.8350535872

FOUND CONC
MEAN= 13.8845 SD= 13.9256762002

NB. RUNS 1 TOTAL X-Y ALL RUNS 20 NB. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.272975323625
SLOPE= 0.735758090615

USE FOR ACCURACY

R= 0.952875566024

MEAN SQR DEV OF POINTS FROM REGRESSION= 18.8379829636

ST ERROR EST= 4.3402745263

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 8.18590696298

x(d)= 21.2428137955

[illegible]

Figure II-40. TNB on Brick - Graph of Target-Found Concentration Points

Table II-43. TNB on Brick - Inaccuracy and Imprecision Data

1,3,5-TRINITROBENZENE (135TNB)

**BRICK SURFACE
STATISTICAL DATA USED TO DETERMINE PERCENT
INACCURACY AND IMPRECISION**

Mn Target Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	1.795	0.543	-28.200	30.257
5.000	4.505	1.424	-9.900	31.608
10.000	7.918	1.309	-20.825	16.534
25.000	17.815	4.692	-28.740	26.335
50.000	37.390	9.231	-25.220	24.689
Means		3.440	-22.577	25.885

1,3,5-TRINITROBENZENE (135TNB)
 BRICK SURFACE
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure II-41. TNB on Brick - Graph of Inaccuracy

1,3,5-TRINITROBENZENE (135TNB)
 BRICK SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-42. TNB on Brick - Graph of Imprecision

Table II-44. TNB on Transite - Target vs. Found Concentrations

1,3,5-TRINITROBENZENE (135TNB)

TRANSITE SURFACE	
TARGET CONC. VS FOUND CONC	
Target Conc ug/10 sq cm	Found Conc ug/10 sq cm
2.500	0.860 2.050 1.780 1.540
5.000	1.600 3.760 4.620 3.730
10.000	3.270 8.230 8.220 8.140
25.000	6.990 17.090 20.220 18.330
50.000	12.180 33.330 34.930 39.400

Table II-45. TNB on Transite - Analysis of Target-
Found Concentration Points

1,3,5-TRINITROBENZENE (135TNB)
TRANSITE SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 18.5 SD= 18.0350535072

FOUND CONC

MEAN= 11.5135 SD= 12.010847849

NB. RUNS 1 TOTAL X-Y ALL RUNS 20 NB. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.355752427104

SLOPE= 0.592310679612

USE FOR ACCURACY

R= 0.889392237864

MEAN SQR DEV OF POINTS FROM REGRESSION= 31.0226365031

ST ERROR EST= 5.64115560706

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 10.840373461

x(d)= 34.6225232479




Table II-46. TNB on Transite - Inaccuracy and Imprecision Data

1,3,5-TRINITROBENZENE (135TNB)

TRANSITE SURFACE
STATISTICAL DATA USED TO DETERMINE PERCENT
INACCURACY AND IMPRECISION

Mn Target Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	1.558	0.510	-37.700	32.715
5.000	3.428	1.286	-31.450	37.529
10.000	6.965	2.464	-30.350	35.372
25.000	15.658	5.920	-37.370	37.889
50.000	29.960	12.128	-40.080	40.482
Means		4.462	-35.390	36.781

1,3,5-TRINITROBENZENE (135TNB)
 TRANSITE SURFACE
 MEAN FOUND CONCENTRATION (REPORT)
 MEAN INACCURACY




Figure II-44. TNB on Transite - Graph of Inaccuracy

1,3,5-TRINITROBENZENE (135TNB)
 TRANSITE SURFACE
 IMPRECISION NEAR FOUND CONCENTRATION (REPORT)




Figure II-45. TNB on Transite - Graph of Imprecision

Table II-47. TNB on Transite (3 days) - Target vs.
Found Concentrations

1,3,5-TRINITROBENZENE (135TNB)

TRANSITE SURFACE	
TARGET CONC. VS. FOUND CONC	
Target Conc ug/10 sq cm	Found Conc ug/10 sq cm
2.500	2.050 1.780 1.540
5.000	3.760 4.620 3.730
10.000	8.230 8.220 8.140
25.000	17.090 20.220 18.330
50.000	33.330 34.930 39.400

Table II-48. TNB on Transite (3 days) - Analysis of Target-Found Concentration Points

1,3,5-TRINITROBENZENE (135TNB)
TRANSITE SURFACE

ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 18.5 SD= 18.1953683274

FOUND CONC

MEAN= 13.691333333 SD= 13.0087014102

N0. RUNS 1 TOTAL X-Y ALL RUNS 15 N0. CONCENTR 15
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.541645091694

SLOPE= 0.710793959008

USE FOR ACCURACY

R= 0.994192846865

MEAN SQR DEV OF POINTS FROM REGRESSION= 2.11040858187

ST ERROR EST= 1.45275207171

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 13

TWO TAIL P LEVEL IS .1

t= 1.77093170942

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

Y(C)= 3.28917417064

X(D)= 7.64791513283

1,3,5-TRINITROBENZENE (135TNB)
TRANSITE SURFACE
FOUND CONC





Figure II-46. TNB on Transite (3 days) - Graph of Target-Found Concentration Points

Table II-49. TNB on Transite (3 days) - Inaccuracy and Imprecision Data

1,3,5-TRINITROBENZENE (135TNB)					
TRANSITE SURFACE					
STATISTICAL DATA USED TO DETERMINE PERCENT					
INACCURACY AND IMPRECISION					
Mn Target Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision	
2.500	1.790	0.255	-28.400	14.254	
5.000	4.037	0.505	-19.267	12.520	
10.000	8.197	0.049	-18.033	0.602	
25.000	18.547	1.576	-25.813	8.499	
50.000	35.887	3.146	-28.227	8.767	
Means		1.106	-23.948	8.928	

1,3,5-TRINITROBENZENE (1351NB)
 TRANSITE SURFACE
 MEAN INACCURACY MEAN FOUND CONCENTRATION <REPORT>



VERTICAL AXIS TIC INTERVAL = 1.03666666667

Figure II-47. TNB on Transite (3 days) - Graph of Inaccuracy

1,3,5-TRINITROBENZENE (135TNB)
 TRANSITE SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-48. TNB on Transite - Graph of Imprecision

Table II-50. 2,4-DNT on Metal - Target vs. Found Concentrations

2,4-DINITROTOLUENE (24DNT)	
METAL SURFACE	
TARGET CONC. us/10 sq cm	US FOUND CONC Found Conc us/10 sq cm
2.500	2.120
	2.300
	2.480
	2.070
5.000	3.710
	4.370
	4.170
	4.840
10.000	9.020
	8.720
	10.040
	11.180
25.000	21.000
	20.820
	25.960
	21.740
50.000	39.160
	38.400
	45.230
	46.510

Table II-51. 2,4-DNT on Metal - Analysis of Target-
Found Concentration Points

2,4-DINITROTOLUENE (24DNT)
METAL SURFACE
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.0350535872

FOUND CONC
MEAN= 16.192 SD= 15.3408483262

NB. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.59741302589
SLOPE= 0.842950647249
USE FOR ACCURACY
R= 0.990992138847
MEAN SQR DEV OF POINTS FROM REGRESSION= 4.45523971932
ST ERROR EST= 2.11974387819
USE FOR PRECISION
T FOR CONFIDENCE BAND
D.F.= 18

TWO TAIL P LEVEL IS .1
t= 1.73406096408
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C
Y(C)= 4.44559634128
X(d)= 9.04498923486

2,4-DINITROTOLUENE (24DNT)
METAL SURFACE
FOUND CONC

46.51




Figure II-49. 2,4-DNT on Metal - Graph of Target-Found Concentration Points

Table II-52. 2,4-DNT on Metal - Inaccuracy and Imprecision Data

2,4-DINITROTOLUENE (24DNT)

METAL SURFACE

STATISTICAL DATA USED TO DETERMINE PERCENT

INACCURACY AND IMPRECISION

Mn Target Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	2.243	0.187	-10.300	8.322
5.000	4.273	0.468	-14.530	10.965
10.000	9.740	1.114	-2.600	11.436
25.000	22.380	2.420	-10.480	10.812
50.000	42.325	4.138	-15.350	9.777
Means		1.665	-10.656	10.263

2,4-DINITROTOLUENE (24DNT)
 METAL SURFACE
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure 11-50. 2,4-DNT on Metal - Graph of Inaccuracy

2,4-DINITROTOLUENE (24DNT)
METAL SURFACE

IMPRECISION MEAN FOUND CONCENTRATION <REPORT>




Figure II-51. 2,4-DNT on Metal - Graph of Imprecision

Table II-53. 2,4-DNT on Concrete - Target vs. Found Concentrations

2,4-DINITROTOLUENE (24DNT) CONCRETE SURFACE TARGET CONC. VS. FOUND CONC.	
Target Conc ug/10 sq cm	Found Conc ug/10 sq cm
2.500	2.530
	2.140
	2.100
	2.090
5.000	4.340
	4.510
	4.140
	3.830
10.000	6.310
	8.360
	7.890
	8.160
25.000	23.120
	20.680
	22.670
	17.900
50.000	33.880
	41.990
	44.240
	36.070

Table II-54. 2,4-DNT on Concrete - Analysis of
Target Found Concentration Points

2,4-DINITROTOLUENE (24DNT)
CONCRETE SURFACE
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.0350535872

FOUND CONC
MEAN= 14.8475 SD= 14.301442281

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.364439724919
SLOPE= 0.782868122977
USE FOR ACCURACY
R= 0.987247878379
MEAN SQR DEV OF POINTS FROM REGRESSION= 5.47110764003
ST ERROR EST= 2.33903989706
USE FOR PRECISION
T FOR CONFIDENCE BAND
D.F.= 18
TWO TAIL P LEVEL IS .1
t= 1.73406096408
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C
Y(C)= 4.62883883032
X(D)= 10.7804175652

Table II-55. 2,4-DNT on Concrete - Inaccuracy and Imprecision Data

2,4-DINITROTOLUENE (24DNT) CONCRETE SURFACE STATISTICAL DATA USED TO DETERMINE PERCENT INACCURACY AND IMPRECISION				
Mn Target Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	2.215	0.211	-11.400	9.531
5.000	4.205	0.292	-15.900	6.948
10.000	7.680	0.933	-23.200	12.154
25.000	21.093	2.378	-15.630	11.273
50.000	39.045	4.871	-21.910	12.476
Means		1.737	-17.608	10.476

2,4-DINITROTOLUENE (24DNT)
 CONCRETE SURFACE
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure II-53. 2,4-DNT on Concrete - Graph of Inaccuracy

2,4-DINITROTOLUENE (24DNT)
 CONCRETE SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION <REPORT>




Figure II-54. 2,4-DNT on Concrete - Graph of Imprecision

Table II-56. 2,4-DNT on Brick - Target vs. Found Concentrations

2,4-DINITROTOLUENE (24DNT) BRICK SURFACE	
TARGET CONC. ug/10 sq cm	US FOUND CONC. Found Conc ug/10 sq cm
2.500	1.530 2.270 1.040 1.780
5.000	2.410 5.420 4.160 4.740
10.000	7.570 6.300 7.420 9.830
25.000	22.860 14.260 13.890 20.880
50.000	37.560 37.440 22.770 42.410

Table II-57. 2,4-DNT on Brick - Analysis of Target-
Found Concentration Points

2,4-DINITROTOLUENE (24DNT)
BRICK SURFACE
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.0350535872

FOUND CONC
MEAN= 13.327 SD= 13.1141246071

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.494765372168
SLOPE= 0.693634304207
USE FOR ACCURACY
R= 0.953912839867
MEAN SQR DEV OF POINTS FROM REGRESSION= 16.347255196
ST ERROR EST= 4.04317390128
USE FOR PRECISION
T FOR CONFIDENCE BAND
D.F.= 18
TWO TAIL P LEVEL IS .1
t= 1.73406096408
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C
Y(C)= 7.86604073574
X(d)= 20.9893901903

Table II-58. 2,4-DNT on Brick - Inaccuracy and Imprecision Data

2,4-DINITROTOLUENE (24DNT) BRICK SURFACE STATISTICAL DATA USED TO DETERMINE PERCENT INACCURACY AND IMPRECISION				
Mn Target Conc ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	1.655	0.512	-33.800	30.961
5.000	4.183	1.289	-16.350	30.819
10.000	7.780	1.479	-22.200	19.016
25.000	17.973	4.575	-28.110	25.455
50.000	35.045	8.505	-29.910	24.267
Means		3.272	-26.074	26.104

2,4-DINITROTOLUENE (24DNT)
 BRICK SURFACE
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure II-56. 2,4-DNT on Brick - Graph of Inaccuracy

2,4-DINITROTOLUENE (24DNT)

BRICK SURFACE

IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-57. 2,4-DNT on Brick - Graph of Imprecision

Table II-59. 2,4-DNT on Transite - Target vs. Found Concentrations

2,4-DINITROTOLUENE (24DNT)	
TRANSITE SURFACE	
TARGET CONC. ug/10 sq cm	US FOUND CONC. Found Conc ug/10 sq cm
2.500	0.700
	2.220
	2.040
	1.660
5.000	1.540
	3.940
	4.250
	3.050
10.000	2.850
	7.620
	7.190
	8.020
25.000	7.950
	16.380
	18.930
	21.220
50.000	11.620
	37.580
	39.370
	41.520

Table II-60. 2,4-DNT on Transite - Analysis of Target-
Found Concentration Points

2,4-DINITROTOLUENE (24DNT)
TRANSITE SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.0350535872

FOUND CONC
MEAN= 11.9825 SD= 13.2377478109

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.0521541262134
SLOPE= 0.65052184466
USE FOR ACCURACY
R= 0.886268306037
MEAN SQR DEV OF POINTS FROM REGRESSION= 39.6820662231
ST ERROR EST= 6.29937030385
USE FOR PRECISION
T FOR CONFIDENCE BAND
D.F.= 18

TWO TAIL P LEVEL IS .1
t= 1.73486096408
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C
(EACH TARGET CONC CONSIDERED INDEP SAMPLE
MEASURED 1 TIME(S))
y(c)= 11.4324849867
x(d)= 35.2284612374




Table II-61. 2,4-DNT on Transite - Inaccuracy and Imprecision Data

2,4-DINITROTOLUENE (24DNT) TRANSITE SURFACE STATISTICAL DATA USED TO DETERMINE PERCENT INACCURACY AND IMPRECISION				
Mn Target Con	Mn Found Conc	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	1.655	0.678	-33.800	48.973
5.000	3.195	1.215	-36.100	38.026
10.000	6.400	2.404	-35.800	37.446
25.000	16.120	5.794	-35.520	35.945
50.000	32.523	14.028	-34.955	43.133
Means		4.824	-35.235	39.104

2,4-DINITROTOLUENE (24DNT)
 TRANSITE SURFACE
 MEAN INACCURACY

MEAN FOUND CONCENTRATION (REPORT)

-33.800

0

-36.10

VERTICAL AXIS TIC INTERVAL = 0.23 TARGET CONC

8

9

10

20

30

40

50

Figure II-59. 2,4-DNT on Transite - Graph of Inaccuracy

2,4-DINITROTOLUENE (24DNT)
 TRANSITE SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-60. 2,4-DNT on Transite - Graph of Imprecision

Table II-62. 2,4-DNT on Transite (3 days) - Target vs.
Found Concentrations

2,4-DINITROTOLUENE (24DNT) TRANSITE SURFACE TARGET CONC. VS. FOUND CONC.	
Target Conc ug/10 sq cm	Found Conc ug/10 sq cm
2.500	2.220 2.040 1.660
5.000	3.940 4.250 3.050
10.000	7.620 7.190 8.020
25.000	16.380 18.930 21.220
50.000	37.580 39.370 41.520

Table 11-63. 2,4-DNT on Transite (3 days) - Analysis
of Target-Found Concentration Points

2,4-DINITROTOLUENE (24DNT)
TRANSITE SURFACE
ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.1953683274

FOUND CONC
MEAN= 14.332666667 SD= 14.4185805061

NØ. RUNS 1 TOTAL X-Y ALL RUNS 15 NØ. CONCENTR 15
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.270698489752
SLOPE= 0.789371089536
USE FOR ACCURACY
R= 0.996138122967
MEAN SQR DEV OF POINTS FROM REGRESSION= 1.72591231227
ST ERROR EST= 1.31373981909
USE FOR PRECISION
T FOR CONFIDENCE BAND
D.F.= 13
TWO TAIL P LEVEL IS .1
t= 1.77893170942
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C
(EACH TARGET CONC CONSIDERED INDEP SAMPLE
MEASURED 1 TIME(S))
Y(C)= 2.21392254463
X(D)= 6.23754930521

2,4-DINITROTOLUENE (24DNT)
TRANSITE SURFACE
FOUND CONC




Figure 11-61. 2,4-DNT on Transite (3 days) - Graph of Target-found Concentration Points.

Table 11-64. 2,4-DNT on Transite (3 days) - Inaccuracy and Imprecision Data

2,4-DINITROTOLUENE (24DNT)
TRANSITE SURFACE

STATISTICAL DATA USED TO DETERMINE PERCENT
INACCURACY AND IMPRECISION

Mn Target Con- ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	1.973	0.286	-21.067	14.488
5.000	3.747	0.623	-25.067	16.626
10.000	7.610	0.415	-23.900	5.455
25.000	18.843	2.421	-24.627	12.849
50.000	39.490	1.973	-21.020	4.996
Means		1.144	-23.136	10.883

2,4-DINITROTOLUENE (24DNT)
 TRANSITE SURFACE
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure 11-62. 2,4-DNT on Transite (3 days) - Graph of Inaccuracy

2,4-DINITROTOLUENE (24DNT)
 TRANSITE SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-63. 2,4-DNT on Transite (3 days) - Graph of Imprecision

Table 11-65. TNT on Metal - Target vs. Found Concentrations

2,4,6-TRINITROTOLUENE (246TNT)

METAL SURFACE	
TARGET CONC. ug/10 sq cm	US FOUND CONC. Found Conc ug/10 sq cm
2.500	1.920
	2.140
	2.120
5.000	1.810
	3.360
	4.510
10.000	4.020
	3.910
25.000	9.560
	9.090
	8.070
50.000	8.680
	23.630
	24.690
	21.180
	20.570
	48.040
	50.620
	43.920
	50.080

2,4,6-TRINITROTOLUENE (246TNT)
METAL SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.0350535872

FOUND CONC
MEAN= 17.096 SD= 17.6143310075

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.90170105178
SLOPE= 0.972848705502
USE FOR ACCURACY
R= 0.996085433419
MEAN SQR DEV OF POINTS FROM REGRESSION= 2.55903491889
ST ERROR EST= 1.59969838373
USE FOR PRECISION
T FOR CONFIDENCE BAND
D.F.= 18

TWO TAIL P LEVEL IS .1
t= 1.73406096408
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C
(EACH TARGET CONC CONSIDERED INDEP SAMPLE
MEASURED 1 TIME(S))
Y(C)= 2.01477436407
X(D)= 5.95490116603

2,4,6-TRINITROTOLUENE (246TNT)
METAL SURFACE
FOUND CONC




Figure 11-64. TNT on Metal - Graph of Target-Found Concentration Points

Table 11-67. TNT on Metal - Inaccuracy and Imprecision Data

2,4,6-TRINITROTOLUENE (246TNT)

METAL SURFACE

STATISTICAL DATA USED TO DETERMINE PERCENT

INACCURACY AND IMPRECISION

Mn Target Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	1.998	0.160	-20.100	7.993
5.000	3.950	0.472	-21.000	11.948
10.000	8.850	0.632	-11.500	7.143
25.000	22.518	1.961	-9.930	8.710
50.000	48.165	3.040	-3.670	6.312
Means		1.253	-13.240	8.421

2,4,6-TRINITROTOLUENE (246TNT)
 METAL SURFACE
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure 11-65. TNT on Metal - Graph of Inaccuracy

2,4,6-TRINITROTOLUENE (246TNT)
 METAL SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-66. TNT on Metal - Graph of Imprecision

Table II-68. TNT on Concrete - Target vs. Found Concentrations

2,4,6-TRINITROTOLUENE (246TNT)

CONCRETE SURFACE	
TARGET CONC. ug/10 sq cm	VS FOUND CONC Found Conc ug/10 sq cm
2.500	1.540
	2.020
	1.310
	1.790
5.000	3.950
	4.430
	3.090
	3.200
10.000	6.960
	7.630
	5.250
	6.690
25.000	19.840
	20.130
	14.290
	17.430
50.000	40.730
	43.930
	27.160
	35.640

Table II-69. TNT on Concrete Analysis of Target-
Found Concentration Points

2,4,6-TRINITROTOLUENE (246TNT)
CONCRETE SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 18.5 SD= 18.0350535872

FOUND CONC

MEAN= 13.3505 SD= 13.7346130078

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.36244538835

SLOPE= 0.741240291262

USE FOR ACCURACY

R= 0.973329817619

MEAN SQR DEV OF POINTS FROM REGRESSION= 10.4794771204

ST ERROR EST= 3.23720205122

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 5.53942979924

x(d)= 15.7269052076

2,4,6-TRINITROTOLUENE (246TNT)
 CONCRETE SURFACE
 FOUND CONC




Figure 11-67. TNT on Concrete - Graph of Target-Found Concentration Points

Table II-70. TNT on Concrete - Inaccuracy and Imprecision Data

2,4,6-TRINITROTOLUENE (246TNT)

CONCRETE SURFACE
STATISTICAL DATA USED TO DETERMINE PERCENT
INACCURACY AND IMPRECISION

Mn Target Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	1.665	0.307	-33.400	18.456
5.000	3.668	0.636	-26.650	17.340
10.000	6.633	1.003	-33.675	15.120
25.000	17.923	2.707	-28.310	15.105
50.000	36.865	7.315	-26.270	19.843
Means		2.394	-29.661	17.173

2,4,6-TRINITROTOLUENE (246TNT)
 CONCRETE SURFACE
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)
 -26.27




Figure 11-68. TNT on Concrete - Graph of Inaccuracy

2,4,6-TRINITROTOLUENE (246TNT)
 CONCRETE SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-69. TNT on Concrete - Graph of Imprecision

Table 11-71. TNT on Brick - Target vs. Found Concentrations

2,4,6-TRINITROTOLUENE (246TNT)

BRICK SURFACE

TARGET CONC. VS. FOUND CONC.

Target Conc ug/10 sq cm	Found Conc ug/10 sq cm
2.500	1.470 1.960 0.890 1.800
5.000	2.320 5.930 4.100 4.180
10.000	8.080 6.280 7.310 9.400
25.000	20.760 15.210 13.750 20.760
50.000	36.650 34.530 24.570 38.460

Table 11-72. TNT on Brick - Analysis of Target-
Found Concentration Points

2,4,6-TRINITROTOLUENE (246TNT)
BRICK SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 18.5 SD= 18.0350535872

FOUND CONC

MEAN= 12.9205 SD= 12.3402646372

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.645345873786

SLOPE= 0.66352184466

USE FOR ACCURACY

R= 0.969724100462

MEAN SQR DEV OF POINTS FROM REGRESSION= 9.50589122507

ST ERROR EST= 3.09610904606

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

Y(C)= 6.28998864488

X(D)= 16.80061912189

2,4,6-TRINITROTOLUENE (246TNT)
BRICK SURFACE
FOUND CONC




Figure 11-70. TNT on Brick - Graph of Target-Found Concentration Points

Table 11-73. TNT on Brick - Inaccuracy and Imprecision Data

2,4,6-TRINITROTOLUENE (246TNT)

BRICK SURFACE STATISTICAL DATA USED TO DETERMINE PERCENT INACCURACY AND IMPRECISION

Mn Target Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	1.530	0.473	-38.800	38.911
5.000	4.133	1.474	-17.350	35.672
10.000	7.768	1.315	-22.325	16.925
25.000	17.620	3.674	-29.520	20.854
50.000	33.553	6.200	-32.895	18.478
Means		2.627	-28.178	24.568

2,4,6-TRINITROTOLUENE (246TNT)
 BRICK SURFACE
 MEAN FOUND CONCENTRATION (REPORT)
 MEAN INACCURACY




Figure 11-71. TNT on Brick - Graph of Inaccuracy

2,4,6-TRINITROTOLUENE (246TNT)
 BRICK SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION <REPORT>




Figure 11-72. TNT on Brick - Graph of Imprecision

Table U-74. TNT on Transite - Target vs. Found Concentrations

2,4,6-TRINITROTOLUENE (246TNT)

TRANSITE SURFACE	
TARGET CONC. ug/10 sq cm	US FOUND CONC Found Conc ug/10 sq cm
2.500	1.160
	1.770
	1.770
	0.900
5.000	1.560
	3.900
	3.510
	2.520
10.000	2.820
	7.770
	6.400
	5.880
25.000	6.450
	17.370
	16.330
	15.800
50.000	11.720
	28.790
	29.470
	32.170

Table 11-75. TNT on Transit - Analysis of Target-Found Concentration Points

2,4,6-TRINITROTOLUENE (246TNT)
TRANSITE SURFACE
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.0350535872

FOUND CONC
MEAN= 9.903 SD= 10.1474591583

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.496468446602

SLOPE= 0.508461165049

USE FOR ACCURACY

R= 0.903686747162

MEAN SQR DEV OF POINTS FROM REGRESSION= 19.9286214266
ST ERROR EST= 4.46414845482

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 8.63523970102

x(d)= 31.8194111464

Arthur D. Little, Inc.

Table 11-76. TNT on Transite - Inaccuracy and Imprecision Data


2,4,6-TRINITROTOLUENE (246TNT)
 TRANSITE SURFACE
 STATISTICAL DATA USED TO DETERMINE PERCENT
 INACCURACY AND IMPRECISION

Mn Targt Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	1.400	0.440	-44.000	31.445
5.000	2.873	1.050	-42.550	36.562
10.000	5.718	2.090	-42.825	36.549
25.000	13.988	5.067	-44.050	36.226
50.000	25.538	9.327	-48.925	36.521
Means		3.595	-44.470	35.461

2,4,6-TRINITROTOLUENE (246TNT)
TRANSITE SURFACE

MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)

-42.55



VERTICAL AXIS TIC INTERVAL = 0.6375

Figure 11-74. TNT on Transite - Graph of Inaccuracy

2,4,6-TRINITROTOLUENE (246TNT)
 TRANSITE SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-75. TNT on Transite - Graph of Imprecision

Table II-77. TNT on Transite (3 days) - Target vs.
Found Concentrations

2,4,6-TRINITROTOLUENE (246TNT)

TRANSITE SURFACE

TARGET CONC. VS. FOUND CONC.	
Target Conc ug/10 sq cm	Found Conc ug/10 sq cm
2.500	1.770
	1.770
	0.900
5.000	3.900
	3.510
	2.520
10.000	7.770
	6.400
	5.880
25.000	17.370
	16.330
	15.800
50.000	28.790
	29.470
	32.170

Table 11-78. TNT on Transite (3 days) - Analysis of
Target-Found Concentration Points

2,4,6-TRINITROTOLUENE (246TNT)

TRANSITE SURFACE

ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 18.5 SD= 18.1953683274

FOUND CONC

MEAN= 11.623333333 SD= 11.0214567787

N0. RUNS 1 TOTAL X-Y ALL RUNS 15 N0. CONCENTR 15
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.468352211435

SLOPE= 0.602971952535

USE FOR ACCURACY

R= 0.995448876481

MEAN SQR DEV OF POINTS FROM REGRESSION= 1.18801497594

ST ERROR EST= 1.08996099744

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 13

TWO TAIL P LEVEL IS .1

t= 1.77093170942

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

Y(C)= 2.52974974353

X(D)= 6.77062537754

32.17

0.00

0 8 0

0 10 20 30 40 50

TARGET CONC

VERTICAL AXIS TIC INTERVAL = 3.217

$y(c)$ $x(d)$

Arthur D Little, Inc

Table 11-79. TNT on Transite (3 days) - Inaccuracy and Imprecision Data

2,4,6-TRINITROTOLUENE (246TNT)

TRANSITE SURFACE

STATISTICAL DATA USED TO DETERMINE PERCENT

INACCURACY AND IMPRECISION

Mn Target Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	1.480	0.502	-40.800	33.939
5.000	3.310	0.711	-33.800	21.493
10.000	6.683	0.976	-33.167	14.689
25.000	16.500	0.799	-34.000	4.841
50.000	30.143	1.788	-39.713	5.931
Means		0.955	-36.296	16.162

2,4,6-TRINITROTOLUENE (246TNT)
 TRANSITE SURFACE
 MEAN INACCURACY




Figure 11-77. TNT on Transite (3 days) - Graph of Inaccuracy

2,4,6-TRINITROTOLUENE (246TNT)
 TRANSITE SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION <REPORT>




Figure 11-78. TNT on Transite (3 days) - Graph of Imprecision

Table 11-80. Tetryl on Metal - Target vs. Found Concentrations

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)

METAL SURFACE	
TARGET CONC. ug/10 sq cm	US FOUND CONC. Found Conc ug/10 sq cm
2.500	1.820
	0.810
	0.000
	0.000
5.000	2.660
	2.430
	0.110
	0.250
10.000	8.520
	8.140
	2.620
	2.070
25.000	22.940
	19.600
	8.610
	11.160
50.000	48.900
	46.740
	41.180
	32.350

Table 11-81. Tetryl on Metal - Analysis of Target-
Found Concentration Points

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)
METAL SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 18.5 SD= 18.0330535872

FOUND CONC

MEAN= 13.0455 SD= 16.5283732046

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -3.22953074434

SLOPE= 0.879731391586

USE FOR ACCURACY

R= 0.959925250553

MEAN SQR DEV OF POINTS FROM REGRESSION= 22.6491360617

ST ERROR EST= 4.7591108478

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C
(EACH TARGET CONC CONSIDERED INDEP SAMPLE

MEASURED 1 TIME(S))

y(c)= 5.44699848777

x(d)= 19.476685492

48.90

0 10 20 30 40 50

0.00

VERTICAL AXIS TIC INTERVAL = 4.89

TARGET CONC

$y(c)$

$x(d)$

Concentration Points

Table II-82. Tetryl on Metal - Inaccuracy and Imprecision Data

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)

METAL SURFACE

STATISTICAL DATA USED TO DETERMINE PERCENT

INACCURACY AND IMPRECISION

Mn Target Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	0.658	0.864	-73.700	131.401
5.000	1.363	1.370	-72.750	100.340
10.000	5.338	3.466	-46.625	64.941
25.000	15.578	6.793	-37.690	43.610
50.000	42.293	7.383	-15.415	17.457
Means		3.975	-49.236	71.590

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)
 METAL SURFACE
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)

-15.42




Figure 11-80. Tetryl on Metal - Graph of Inaccuracy

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)
 METAL SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-81. Tetryl on Metal - Graph of Imprecision

Table II-83. Tetryl on Concrete - Target vs. Found Concentrations

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)

TETRYL CONC. VS. FOUND CONC.	
Target Conc ug/10 sq cm	Found Conc ug/10 sq cm
2.500	0.230
	0.520
	0.250
	0.350
5.000	1.500
	1.790
	0.350
	1.100
10.000	2.880
	3.960
	1.280
	1.250
25.000	15.930
	14.650
	3.460
	8.500
50.000	28.080
	40.160
	9.780
	18.020

Table 11-84. Tetryl on Concrete - Analysis of
Target-Found Concentration Points

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)
CONCRETE SURFACE
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.0350535872

FOUND CONC
MEAN= 7.702 SD= 10.8150614277

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -1.68730380259

SLOPE= 0.507529935275

USE FOR ACCURACY

R= 0.846350216412

MEAN SQR DEV OF POINTS FROM REGRESSION= 35.0255619146
ST ERROR EST= 5.91823976488

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

(EACH TARGET CONC CONSIDERED INDEP SAMPLE

MEASURED 1 TIME(S))

y(c)= 9.10248066865

x(d)= 42.9096479565

Figure 11-82. Tetraol on Concrete - Graph of Target - Found Concentration Points

Table 11-35. Tetryl on Concrete - Inaccuracy and Imprecision Data

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)

**CONCRETE SURFACE
STATISTICAL DATA USED TO DETERMINE PERCENT
INACCURACY AND IMPRECISION**

Mn Target Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	0.338	0.133	-86.500	39.262
5.000	1.185	0.624	-76.300	52.694
10.000	2.343	1.320	-76.575	56.353
25.000	10.635	5.779	-57.460	54.341
50.000	24.010	13.112	-51.980	54.610
Means		4.194	-69.763	51.452

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)
 CONCRETE SURFACE
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure 11-83. Tetryl on Concrete - Graph of Inaccuracy

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)
 CONCRETE SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-84. Tetryl on Concrete - Graph of Imprecision

Table 11-86. Tetryl on Brick - Target vs. Found Concentrations

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)

TARGET CONC. VS. FOUND CONC.	
Target Conc ug/10 sq cm	Found Conc ug/10 sq cm
2.500	2.210
	1.280
	0.620
	0.060
5.000	2.300
	6.050
	3.910
	0.600
10.000	9.340
	6.520
	8.060
	2.560
25.000	24.170
	17.540
	17.820
	11.820
50.000	41.510
	41.820
	30.920
	20.310

Table 11-87. Tetryl on Brick - Analysis of Target-
Found Concentration Points

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)
BRICK SURFACE
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.0350535872

FOUND CONC
MEAN= 12.471 SD= 13.2601609006

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.212591019417

SLOPE= 0.685599514563

USE FOR ACCURACY

R= 0.931916945925

MEAN SQR DEV OF POINTS FROM REGRESSION= 24.441622711

ST ERROR EST= 4.94384695465

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

(EACH TARGET CONC CONSIDERED INDEP SAMPLE

MEASURED 1 TIME(S))

y(c)= 8.80073814918

x(d)= 26.0154833108

Table 11-88. Tetryl on Brick - Inaccuracy and Imprecision Data

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)

BRICK SURFACE
STATISTICAL DATA USED TO DETERMINE PERCENT
INACCURACY AND IMPRECISION

Mn Targt Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	1.043	0.924	-58.300	88.667
5.000	3.215	2.323	-35.700	72.270
10.000	6.620	2.942	-33.800	44.441
25.000	17.838	5.046	-28.650	28.291
50.000	33.640	10.230	-32.720	30.409
Means		4.293	-37.834	52.816

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)
 BRICK SURFACE
 MEAN INACCURACY




Figure 11-86. Tetryl on Brick - Graph of Inaccuracy

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)
 BRICK SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-S7. Tetryl on Brick - Graph of Imprecision

Table 11-89. Tetryl on Transite - Target vs. Found Concentrations

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)

TRANSITE SURFACE	
TARGET CONC. US	FOUND CONC
ug/10 sq cm	ug/10 sq cm
2.500	0.950 1.340 0.140 0.170
5.000	1.460 3.080 0.540 0.450
10.000	2.670 7.850 2.640 0.900
25.000	6.880 22.300 9.000 4.300
50.000	12.710 40.460 21.230 35.070

Table II-90. Tetryl on Transite - Analysis of
Target-Found Concentration Points

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)
TRANSITE SURFACE
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.0350535872

FOUND CONC
MEAN= 8.711 SD= 11.9174652462

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -1.73395631068
SLOPE= 0.56459223301

USE FOR ACCURACY

R= 0.854414170035

MEAN SQR DEV OF POINTS FROM REGRESSION= 40.4738695793
ST ERROR EST= 6.36190769968

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C
(EACH TARGET CONC CONSIDERED INDEP SAMPLE

MEASURED 1 TIME(S))

Y(C)= 9.86469728103

X(D)= 41.354710617

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)
 TRANSITE SURFACE
 FOUND CONC




Figure 11-88. Tetryl on Transite - Graph of Target-Found Concentration Points

Table 11-91. Tetryl on Transite - Inaccuracy and Imprecision Data

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)
 TRANSITE SURFACE
 STATISTICAL DATA USED TO DETERMINE PERCENT
 INACCURACY AND IMPRECISION

Mn Target Conc ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	0.650	0.593	-74.000	91.302
5.000	1.383	1.220	-72.350	88.262
10.000	3.515	3.006	-64.850	85.522
25.000	10.640	8.015	-57.440	75.329
50.000	27.368	12.692	-45.265	46.376
Means		5.105	-62.781	77.358

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)
 TRANSITE SURFACE
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure 11-89. Tetryl on Transite - Graph of Inaccuracy

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)
 TRANSITE SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-90. Tetryl on Transite - Graph of Imprecision

Table II-92. DPA on Metal - Target vs. Found Concentrations

DIPHENYLAMINE (DPA) METAL SURFACE	
TARGET CONC. Target Conc ng/ml	VS FOUND CONC Found Conc ng/ml
2.500	2.050
	2.610
	2.560
	2.410
5.000	3.770
	5.070
	4.980
	4.720
10.000	9.520
	10.140
	10.270
	9.910
25.000	25.240
	25.420
	24.400
	21.800
50.000	49.370
	47.940
	47.220
	43.510

Table II-93. DPA on Metal - Analysis of Target-Found Concentration Points

DIPHENYLAMINE (DPA)
METAL SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.0350535872

FOUND CONC
MEAN= 17.6455 SD= 17.002404621

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.253584142395
SLOPE= 0.940103559871

USE FOR ACCURACY

R= 0.997201187585

MEAN SQR DEV OF POINTS FROM REGRESSION= 1.70567937344
ST ERROR EST= 1.30601660534

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C
(EACH TARGET CONC CONSIDERED INDEP SAMPLE

MEASURED 1 TIME(S))

y(c)= 2.63463632131

x(d)= 5.03549506469




Table II-94. DPA on Metal - Inaccuracy and Imprecision Data

DIPHENYLAMINE (DPA) METAL SURFACE STATISTICAL DATA USED TO DETERMINE PERCENT INACCURACY AND IMPRECISION				
Mn Target Con ng/ml	Mn Found Conc ng/ml	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	2.400	0.253	-3.700	10.510
5.000	4.635	0.595	-7.300	12.847
10.000	9.950	0.329	-0.400	3.303
25.000	24.215	1.670	-3.140	6.898
50.000	47.010	2.499	-5.980	5.315
Means		1.069	-4.104	7.774




Figure II-92. DPA on Metal - Graph of Inaccuracy

DIPHENYLAMINE (DPA)
METAL SURFACE
IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-93. DPA on Metal - Graph of Imprecision

Table 11-95. DPA on Concrete - Target vs. Found Concentrations

DIPHENYLAMINE (DPA) CONCRETE SURFACE TARGET CONC. VS. FOUND CONC.	
Target Conc ng/ml	Found Conc ng/ml
2.500	2.010
	2.200
	1.720
	1.560
5.000	4.940
	4.130
	4.140
	2.950
10.000	5.300
	8.470
	7.470
	4.880
25.000	19.400
	21.100
	19.580
	9.750
50.000	35.410
	43.400
	47.480
	22.730

Table 11-96. DPA on Concrete - Analysis of Target-
Found Concentration Points

DIPHENYLAMINE (DPA)
CONCRETE SURFACE
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.0350535072

FOUND CONC
MEAN= 13.431 SD= 14.2563485027

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.303319174757
SLOPE= 0.742395631068
USE FOR ACCURACY
R= 0.939170713081
MEAN SQR DEV OF POINTS FROM REGRESSION= 25.3061729268
ST ERROR EST= 5.03052412049
USE FOR PRECISION
T FOR CONFIDENCE BAND
D.F.= 18
TWO TAIL P LEVEL IS .1
t= 1.73406096408
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C
(EACH TARGET CONC CONSIDERED INDEP SAMPLE
MEASURED 1 TIME(S))
Y(C)= 8.86803467264
X(D)= 24.4266262877

DIPHENYLAMINE (DPA)
CONCRETE SURFACE
FOUND CONC




Figure 11-94. DPA on Concrete - Graph of Target - Found Concentration Points

Table II-97. DPA on Concrete - Inaccuracy and Imprecision Data

DIPHENYLAMINE (DPA) CONCRETE SURFACE STATISTICAL DATA USED TO DETERMINE PERCENT INACCURACY AND IMPRECISION				
Mn Target Conc ng/ml	Mn Found Conc ng/ml	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	1.873	0.287	-25.100	15.326
5.000	4.040	0.820	-19.200	20.292
10.000	6.530	1.721	-34.700	26.351
25.000	17.450	5.195	-30.170	29.756
50.000	37.255	10.904	-25.490	29.269
Means		3.785	-26.932	24.199

DIPHENYLAMINE (DPA)
 CONCRETE SURFACE
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)

-19.20

-34.70

VERTICAL AXIS TIC INTERVAL = 1.55 TARGET CONC




Figure 11-95. DPA on Concrete - Graph of Inaccuracy

DIPHENYLAMINE (DPA)
CONCRETE SURFACE
IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-96. DPA on Concrete - Graph of Imprecision

Table II-98. BPA on Concrete (3 days) - Target vs. Found Concentrations

DIPHENYLAMINE (DPA) CONCRETE SURFACE TARGET CONC. VS FOUND CONC	
Target Conc ng/ml	Found Conc ng/ml
2.500	2.010
	2.200
	1.720
5.000	4.940
	4.130
	4.140
10.000	5.300
	8.470
	7.470
25.000	19.400
	21.100
	19.580
50.000	35.410
	43.400
	47.480

Table 11-99. BPA on Concrete (3 days) - Analysis of Target-Found Concentration Points

DIPHENYLAMINE (DPA)
CONCRETE SURFACE

ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.1953683274

FOUND CONC
MEAN= 15.1166666667 SD= 15.5726285992

NO. RUNS 1 TOTAL X-Y ALL RUNS 15 NO. CONCENTR 15
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.509147788565

SLOPE= 0.844638619202

USE FOR ACCURACY

R= 0.986893039028

MEAN SQR DEV OF POINTS FROM REGRESSION= 6.80118484773
ST ERROR EST= 2.60790813637

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F. = 13

TWO TAIL P LEVEL IS .1

t= 1.77093170942

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

(EACH TARGET CONC CONSIDERED INDEP SAMPLE

MEASURED 1 TIME(S))

Y(C)= 4.42307951945

X(D)= 11.514520374




Table 11-100. DPA on Concrete (3 days) - Inaccuracy and Imprecision Data

DIPHENYLAMINE (DPA)

CONCRETE SURFACE

STATISTICAL DATA USED TO DETERMINE PERCENT

INACCURACY AND IMPRECISION

Mn Target Conc ng/ml	Mn Found Conc ng/ml	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	1.977	0.242	-20.933	12.229
5.000	4.403	0.465	-11.933	10.556
10.000	7.080	1.621	-29.200	22.890
25.000	20.027	0.934	-19.893	4.663
50.000	42.097	6.140	-15.807	14.585
Means		1.880	-19.553	12.984

DIPHENYLAMINE (DPA)
 CONCRETE SURFACE
 MEAN INACCURACY
 -11.93
 MEAN FOUND CONCENTRATION (REPORT)




Figure 11-98. DPA on Concrete (3 days) - Graph of Inaccuracy

DIPHENYLAMINE (DPA)
 CONCRETE SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-99. BPA on Concrete (3 days) - Graph of Imprecision

Table II-101. DPA on Brick - Target vs. Found Concentrations

DIPHENYLAMINE (DPA) BRICK SURFACE	
TARGET CONC. ng/ml	VS. FOUND CONC. Found Conc ng/ml
2.500	1.920
	2.180
	1.710
	2.200
5.000	2.260
	4.950
	4.550
	4.750
10.000	7.900
	7.060
	7.890
	8.260
25.000	22.880
	15.150
	13.180
	16.430
50.000	35.990
	44.510
	23.090
	38.060

Table II-102. DPA on Brick - Analysis of Target-
Found Concentration Points

DIPHENYLAMINE (DPA)
BRICK SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.0350535872

FOUND CONC
MEAN= 13.246 SD= 13.1338996494

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.41527710356
SLOPE= 0.693552588997

USE FOR ACCURACY

R= 0.952364373267

MEAN SQR DEV OF POINTS FROM REGRESSION= 16.9340657172
ST ERROR EST= 4.11510215148

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C
(EACH TARGET CONC CONSIDERED INDEP SAMPLE

MEASURED 1 TIME(S))

Y(C)= 7.91768779587

X(D)= 21.366954558

DIPHENYLAMINE (DPA)
BRICK SURFACE
FOUND CONC

44.51




Figure 11-10b. DPA on Brick - Graph of Target-Found Concentration Points

Table 11-103. DPA on Brick - Inaccuracy and Imprecision Data

DIPHENYLAMINE (DPA)

BRICK SURFACE

STATISTICAL DATA USED TO DETERMINE PERCENT INACCURACY AND IMPRECISION

Mn Target Con ng/ml	Mn Found Conc ng/ml	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	2.003	0.233	-19.900	11.636
5.000	4.128	1.256	-17.450	36.422
10.000	7.778	0.508	-22.225	6.536
25.000	16.910	4.198	-32.360	24.828
50.000	35.413	8.981	-29.175	25.360
Means		3.835	-24.222	19.756

DIPHENYLAMINE (DPA)
 BRICK SURFACE
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure 11-101. DPA on Brick - Graph of Inaccuracy

DIPHENYLAMINE (DPA)
BRICK SURFACE
IMPRECISION MEAN FOUND CONCENTRATION (REPORT)

30.42




Figure II-102. DPA on Brick - Graph of Imprecision

Table II-104. DPA on Transite - Target vs. Found Concentrations

DIPHENYLAMINE (DPA) TRANSITE SURFACE TARGET CONC. VS. FOUND CONC.	
Target Conc ng/ml	Found Conc ng/ml
2.500	0.810
	1.490
	1.880
	1.740
5.000	1.320
	4.640
	4.740
	3.910
10.000	2.840
	9.790
	9.030
	7.970
25.000	7.550
	18.970
	21.120
	17.340
50.000	11.650
	45.170
	37.090
	41.580

Table II-105. DPA on Transite - Analysis of Target-
Found Concentration Points

DIPHENYLAMINE (DPA)
TRANSITE SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 18.5 SD= 18.0350535872

FOUND CONC
MEAN= 12.5315 SD= 13.8086106908

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.11885315534

SLOPE= 0.670953883495

USE FOR ACCURACY

R= 0.876314751306

MEAN SQR DEV OF POINTS FROM REGRESSION= 46.7094406587

ST ERROR EST= 6.83443052922

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

(EACH TARGET CONC CONSIDERED INDEP SAMPLE

MEASURED 1 TIME(S))

y(c)= 12.5789824038

x(d)= 37.1490430551

DIPHENYLAMINE (DPA)
TRANSITE SURFACE
FOUND CONC

45.17




Figure 11-103. DPA on Transite - Graph of Target-Found Concentration Points

Table II-106. DPA on Transite - Inaccuracy and Imprecision Data

DIPHENYLAMINE (DPA) TRANSITE SURFACE STATISTICAL DATA USED TO DETERMINE PERCENT INACCURACY AND IMPRECISION				
Mn Target Con ng/ml	Mn Found Conc ng/ml	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	1.480	0.475	-40.800	32.080
5.000	3.653	1.590	-26.950	43.762
10.000	7.408	3.135	-25.925	42.324
25.000	16.245	6.000	-35.020	36.933
50.000	33.873	15.179	-32.255	44.813
Means		5.278	-32.190	39.984

DIPHENYLAMINE (DPA)
 TRANSITE SURFACE
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)

-25.93




Figure 11-104. DPA on Transite - Graph of Inaccuracy

DIPHENYLAMINE (DPA)
 TRANSITE SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-105. DPA on Transite - Graph of Imprecision

Table II-107. DPA on Transite (3 days) - Target vs. Found Concentrations

DIPHENYLAMINE (DPA) TRANSITE SURFACE TARGET CONC. VS. FOUND CONC.	
Target Conc ng/ml	Found Conc ng/ml
2.500	1.490
	1.880
	1.740
5.000	4.640
	4.740
	3.910
10.000	9.790
	9.030
	7.970
25.000	18.970
	21.120
	17.340
50.000	45.170
	37.090
	41.580

Table II-108. DPA on Transite (3 days) - Target vs.
Found Concentrations

DIPHENYLAMINE (DPA)

TRANSITE SURFACE

ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 18.5 SD= 18.1953683274

FOUND CONC

MEAN= 15.097333333 SD= 14.9797851723

N0. RUNS 1 TOTAL X-Y ALL RUNS 15 N0. CONCENTR 15

MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.011459007551

SLOPE= 0.816691477886

USE FOR ACCURACY

R= 0.992803695584

MEAN SQR DEV OF POINTS FROM REGRESSION= 3.84924285745

ST ERROR EST= 1.96194873976

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 13

TWO TAIL P LEVEL IS .1

t= 1.77093178942

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

(EACH TARGET CONC CONSIDERED INDEP SAMPLE

MEASURED 1 TIME(S))

y(c)= 3.69909237891

x(d)= 8.97736979202

Table II-109. DPA on Transite (3 days) - Inaccuracy and Imprecision Data


DIPHENYLAMINE (DPA)

**TRANSITE SURFACE
STATISTICAL DATA USED TO DETERMINE PERCENT
INACCURACY AND IMPRECISION**

Mn Target Con ng/ml	Mn Found Conc ng/ml	Standard Deviation	Mean Pct Inaccuracy	Imprecision
2.500	1.703	0.198	-31.867	11.599
5.000	4.430	0.453	-11.400	10.228
10.000	8.930	0.914	-10.700	10.236
25.000	19.143	1.896	-23.427	9.984
50.000	41.280	4.848	-17.440	9.807
Means		1.502	-18.967	10.355

DIPHENYLAMINE (DPA)
 TRANSITE SURFACE
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)

-10.70



-31.87

TARGET CONC

VERTICAL AXIS TIC INTERVAL = 2.11666666667

Figure 11-107. DPA on Transite (3 days) - Graph of Inaccuracy




Figure 11-108. DPA on Transite (3 days) - Graph of Imprecision

Table II-110. 2,6-DNT on Metal - Target vs. Found Concentrations

2,6-DINITROTOLUENE (26DNT) METAL SURFACE	
TARGET CONC. ug/10 sq cm	VS. FOUND CONC. Found Conc ug/10 sq cm
10.000	9.550 7.500 11.350 10.500
20.000	19.000 21.000 22.500 19.300
40.000	29.000 39.900 40.000 37.100
100.000	86.700 97.000 95.900 93.600
200.000	168.000 208.000 192.000 189.000

Table II-111 - 2,6-DNT on Metal - Analysis of Target-
Found Concentration Points

2,6-DINITROTOLUENE (26DNT)
METAL SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 74 SD= 72.140214349

FOUND CONC

MEAN= 69.845 SD= 68.3900825739

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.0850323624595

SLOPE= 0.942702265372

USE FOR ACCURACY

R= 0.994394814712

MEAN SQR DEV OF POINTS FROM REGRESSION= 55.1910252609
ST ERROR EST= 7.42906624421

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 13.6292662533

x(d)= 28.5092514656

2,6-DINITROTOLUENE (26DNT)
METAL SURFACE
FOUND CONC




Figure II-109. 2,6-DNT on Metal - Graph of Target-
Found Concentration Points

TABLE II-112. 2,6-DNT on Metal - Inaccuracy and Imprecision Data

2,6-DINITROTOLUENE (26DNT)
METAL SURFACE
STATISTICAL DATA USED TO DETERMINE PERCENT
INACCURACY AND IMPRECISION

Mn Target Conc ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
10.000	9.725	1.656	-2.750	17.024
20.000	20.450	1.626	2.250	7.950
40.000	36.500	5.178	-8.750	14.185
100.000	93.300	4.622	-6.700	4.954
200.000	189.250	16.439	-5.375	11.687
Means		5.904	-4.265	10.560

2,6-DINITROTOLUENE (26DNT)
 METAL SURFACE
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure 11-110. 2,6-DNT on Metal - Graph of Inaccuracy

2,6-DINITROTOLUENE (26DNT)
 METAL SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-111. 2,6-DNT on Metal - Graph of Imprecision

Table II-113. 2,6-DNT on Concrete - Target vs.
Found Concentrations

2,6-DINITROTOLUENE (26DNT)	
CONCRETE	
TARGET CONC. US FOUND CONC	
Target Conc ug/10 sq cm	Found Conc ug/10 sq cm
10.000	15.000
	10.200
	9.000
	10.200
20.000	10.000
	16.600
	19.000
	20.400
40.000	37.600
	35.000
	37.400
	35.700
100.000	76.900
	94.600
	97.200
	104.000
200.000	120.000
	189.000
	181.000
	196.000

Table II-114. 2,6-DNT on Concrete - Analysis of Target-Found Concentration Points

2,6-DINITROTOLUENE (26DNT)
CONCRETE SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 74 SD= 72.140214349

FOUND CONC
MEAN= 65.53 SD= 63.8466262535

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 1.8729368932

SLOPE= 0.860230582524

USE FOR ACCURACY

R= 0.971973340682

MEAN SQR DEV OF POINTS FROM REGRESSION= 237.809596818

ST ERROR EST= 15.4210763832

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 29.9877306664

x(d)= 64.5501819952

Table II-115. 2,6-DNT on Concrete - Inaccuracy and Imprecision Data

2,6-DINITROTOLUENE (26DNT)

CONCRETE STATISTICAL DATA USED TO DETERMINE PERCENT INACCURACY AND IMPRECISION

Mn Target Conc ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
10.000	9.850	0.574	-1.500	5.832
20.000	16.500	4.609	-17.500	27.931
40.000	36.625	1.014	-8.438	2.770
100.000	93.175	11.551	-6.825	12.397
200.000	171.500	34.876	-14.250	20.336
Means		10.525	-9.702	13.853

2,6-DINITROTOLUENE (26DNT)
 CONCRETE
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure II-113. 2,6-DNT on Concrete - Graph of Inaccuracy

2,6-DINITROTOLUENE (26DNT)
CONCRETE
IMPRECISION MEAN FOUND CONCENTRATION (REPORT)

27.93




Figure 11-114. 2,6-DNT on Concrete - Graph of Imprecision

Table II-114. 2,6-DNT on Brick - Target vs.
Found Concentrations

2,6-DINITROTOLUENE (26DNT) BRICK SURFACE	
TARGET CONC. ug/10 sq cm	VS FOUND CONC. Found Conc ug/10 sq cm
10.000	4.000 5.300 5.600 4.300
20.000	8.000 9.100 8.680 7.100
40.000	17.300 24.800 25.200 16.700
100.000	46.200 62.400 57.400 45.400
200.000	89.800 94.300 123.000 80.820

TABLE II-117. 2,6-DNT on Brick - Analysis of Target-
Found Concentration Points

2,6-DINITROTOLUENE (26DNT)
BRICK SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 74 SD= 72.140214349

FOUND CONC
MEAN= 36.77 SD= 36.3865408894

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.463001618123

SLOPE= 0.490635113269

USE FOR ACCURACY

R= 0.972736659578

MEAN SQR DEV OF POINTS FROM REGRESSION= 75.1641619382
ST ERROR EST= 8.66972675107

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 16.2691343172

x(d)= 63.6301961361

2,6-DINITROTOLUENE (26DNT)
BRICK SURFACE
FOUND CONC




Figure 11-115. 2,6-DNT on Brick - Graph of Target-Found Concentration Points

TABLE 11-118. 2,6-DNT on Brick - Inaccuracy and Imprecision Data

2,6-DINITROTOLUENE (26DNT)
BRICK SURFACE
STATISTICAL DATA USED TO DETERMINE PERCENT
INACCURACY AND IMPRECISION

Mn Target Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
10.000	4.800	0.770	-52.000	16.040
20.000	0.220	0.073	-50.900	10.626
40.000	21.000	4.620	-47.500	22.039
100.000	52.850	8.399	-47.150	15.092
200.000	96.980	10.229	-51.510	10.797
Means		6.580	-51.412	16.680

2,6-DINITROTOLUENE (26DNT)
BRICK SURFACE
MEAN INACCURACY

-47.15

MEAN FOUND CONCENTRATION (REPORT)




Figure 11-146. 2,6-DNT on Brick - Graph of Inaccuracy

2,6-DINITROTOLUENE (26DNT)
BRICK SURFACE
IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-117. 2,6-DNT on Brick - Graph of Imprecision

TABLE II-119. 2,6-DNT on Transite - Target vs.
Found Concentrations

2,6-DINITROTOLUENE (26DNT)	
TRANSITE SURFACE	
TARGET CONC. VS. FOUND CONC.	
Target Conc ug/10 sq cm	Found Conc ug/10 sq cm
10.000	7.500
	8.700
	7.500
	7.800
20.000	14.700
	16.300
	17.800
	19.900
40.000	27.700
	36.000
	36.600
	33.200
100.000	89.000
	87.300
	81.900
	89.700
200.000	152.000
	156.000
	164.000
	160.000

TABLE II-120. 2,6-DNT on Transite - Analysis of
Target-Found Concentration Points

2,6-DINITROTOLUENE (26DNT)
TRANSITE SURFACE
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 74 SD= 72.140214349

FOUND CONC
MEAN= 60.68 SD= 57.3720507785

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 1.99743122977

SLOPE= 0.793007686084

USE FOR ACCURACY

R= 0.99713612601

MEAN SQR DEV OF POINTS FROM REGRESSION= 19.0720941201
ST ERROR EST= 4.45781158418

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

Y(C)= 10.1246494767

X(D)= 20.3745583959

2,6-DINITROTOLUENE (26DNT)
TRANSITE SURFACE
FOUND CONC





Figure 11-118. 2,6-DNT on Transite - Graph of Target-Found Concentration Points

TABLE II-121. 2,6-DNT on Transite - Inaccuracy and Imprecision Data

2,6-DINITROTOLUENE (26DNT)
 TRANSITE SURFACE
 STATISTICAL DATA USED TO DETERMINE PERCENT
 INACCURACY AND IMPRECISION

Mn Target Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
10.000	7.875	0.560	-21.250	7.211
20.000	17.175	2.214	-14.125	12.892
40.000	33.375	4.063	-16.563	12.174
100.000	86.975	3.530	-13.025	4.059
200.000	158.000	5.164	-21.000	3.260
Means		3.108	-17.193	7.921

2,6-DINITROTOLUENE (26DNT)
 TRANSITE SURFACE
 MEAN FOUND CONCENTRATION (REPORT)
 MEAN INACCURACY
 -13.03



-21.25
 VERTICAL AXIS TIC INTERVAL = 0.8225

Figure II-119. 2,6-DNT on Transite - Graph of Inaccuracy

2,6-DINITROTOLUENE (26DNT)
 TRANSITE SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-120. 2,6-DNT on Transite - Graph of Imprecision

TABLE II-122. NG on Metal - Target vs. Found Concentrations

NITROGLYCERIN METAL SURFACE TARGET CONC. VS. FOUND CONC.	
Target Conc ug/10 sq cm	Found Conc ug/10 sq cm
125.000	122.000 99.500 102.000 119.000
250.000	205.000 210.000 206.000 220.000
500.000	415.000 423.000 402.000 443.000
1250.000	1150.000 1200.000 1150.000 1140.000

TABLE II-123. NG on Metal - Analysis of Target-
Found Concentration Points

NITROGLYCERINE (NG)
METAL SURFACE

ANALYSIS OF 16 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 531.25 SD= 450.693909433

FOUND CONC

MEAN= 475.40625 SD= 424.537778308

NO. RUNS 1 TOTAL X-Y ALL RUNS 16 NO. CONCENTR 16
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -24.3256419256

SLOPE= 0.940671794872

USE FOR ACCURACY

R= 0.998580335082

MEAN SQR DEV OF POINTS FROM REGRESSION= 547.954235351
ST ERROR EST= 23.4084223166

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 14

TWO TAIL P LEVEL IS .1

t= 1.7613101065

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 19.9865337733

x(d)= 93.6030226134




Figure 11-121. NG on Metal - Graph of Target-Found Concentration Points

TABLE II-124. NG on Metal - Inaccuracy and Imprecision Data

NITROGLYCERIN (NG) METAL SURFACE STATISTICAL DATA USED TO DETERMINE PERCENT INACCURACY AND IMPRECISION				
Mn Target Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
125.000	110.625	11.514	-11.500	10.408
250.000	210.250	6.050	-15.900	3.258
500.000	420.750	17.173	-15.050	4.082
1250.000	1160.000	27.000	-7.200	2.334
Means		15.654	-12.613	5.020

NITROGLYCERIN (NG)
 METAL SURFACE
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure 11-122. NG on Metal - Graph of Inaccuracy

NITROGLYCERIN (NG)
METAL SURFACE

IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-123. NG on Metal - Graph of Imprecision

TABLE II-125. NG on Concrete Target vs. Found Concentrations

NITROGLYCERIN (NG) CONCRETE SURFACE TARGET CONC. VS FOUND CONC.	
Target Conc ug/10 sq cm	Found Conc ug/10 sq cm
125.000	93.300
	86.600
	95.900
	99.500
250.000	181.000
	142.000
	193.000
	196.000
500.000	396.000
	392.000
	436.000
	429.000
1250.000	866.000
	883.000
	1007.000
	1014.000

TABLE II-126. NG on Concrete - Analysis of Target-
Found Concentration Points

NITROGLYCERINE (NG)

CONCRETE SURFACE

ANALYSIS OF 16 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 531.25 SD= 450.693909433

FOUND CONC

MEAN= 406.89375 SD= 343.646924364

NO. RUNS 1 TOTAL X-Y ALL RUNS 16 NO. CONCENTR 16
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 4.94128205129

SLOPE= 0.756616410256

USE FOR ACCURACY

R= 0.992304553607

MEAN SQR DEV OF POINTS FROM REGRESSION= 1939.89261997

ST ERROR EST= 44.0442121052

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 14

TWO TAIL P LEVEL IS .1

t= 1.7613101065

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

Y(C)= 88.3170343321

X(C)= 217.473873833

NITROGLYCERINE (NG)
CONCRETE SURFACE
FOUND CONC

1014.00




Figure II-124. NG on Concrete - Graph of Target-Found Concentration Points

TABLE II-127. NG on Concrete - Inaccuracy and Imprecision Data

**NITROGLYCERIN (NG)
CONCRETE SURFACE
STATISTICAL DATA USED TO DETERMINE PERCENT
INACCURACY AND IMPRECISION**

Mn Target Conc ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Fct Inaccuracy	Imprecision
125.000	93.825	5.446	-24.940	5.885
250.000	178.000	24.860	-28.800	13.966
500.000	413.250	22.470	-17.350	5.437
1250.000	942.500	78.878	-24.600	8.369
Means		32.913	-23.923	8.394

NITROGLYCERIN (NG)
 CONCRETE SURFACE
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure II-125. NG on Concrete - Graph of Inaccuracy

NITROGLYCERIN (NG)
 CONCRETE SURFACE
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-126. NG on Concrete - Graph of Imprecision

TABLE II-128. NG on Brick - Target vs. Found Concentrations

NITROGLYCERIN BRICK SURFACE TARGET CONC. VS. FOUND CONC.	
Target Conc ug/10 sq cm	Found Conc ug/10 sq cm
125.000	33.800
	56.600
	58.800
	42.800
250.000	85.900
	87.500
	88.800
	77.400
500.000	222.000
	282.000
	242.000
	202.000
1250.000	537.000
	614.000
	560.000
	453.000

TABLE II-129. NG on Brick - Analysis of Target-
Found Concentration Points

NITROGLYCERINE (NG)
BRICK SURFACE
ANALYSIS OF 16 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 531.25 SD= 450.693909433

FOUND CONC
MEAN= 227.725 SD= 203.495737875

NO. RUNS 1 TOTAL X-Y ALL RUNS 16 NO. CONCENTR 16
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -8.04974358974
SLOPE= 0.443811282051
USE FOR ACCURACY
R= 0.982934796801
MEAN SQR DEV OF POINTS FROM REGRESSION= 1501.39082051
ST ERROR EST= 38.7477847175
USE FOR PRECISION
T FOR CONFIDENCE BAND
D.F. = 14
TWO TAIL P LEVEL IS .1
t= 1.7613101065
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C
Y(C)= 65.29986478
X(D)= 324.819610735

NITROGLYCERINE (NG)
BRICK SURFACE
FOUND CONC

614.00




Figure II-127. NG on Brick - Graph of Target-Pound Concentration Points

TABLE II-130. NG on Brick - Inaccuracy and Imprecision Data

NITROGLYCERIN (NG)					
BRICK SURFACE					
STATISTICAL DATA USED TO DETERMINE PERCENT					
INACCURACY AND IMPRECISION					
Mn Target Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision	
125.000	40.000	11.822	-61.600	24.629	
250.000	84.900	5.139	-66.040	6.053	
500.000	237.000	34.157	-52.600	14.412	
1250.000	541.000	66.958	-56.720	12.377	
Means		29.519	-59.240	14.368	




Figure 11-128. NG on Brick - Graph of Inaccuracy




Figure II-129. NC on Brick - Graph of Imprecision

TABLE II-131. NC on Transite - Target vs. Found Concentrations

NITROGLYCERIN (NG) TRANSITE SURFACE TARGET CONC. VS FOUND CONC.	
Target Conc ug/10 sq cm	Found Conc ug/10 sq cm
125.000	83.600
	76.000
	91.650
	85.500
250.000	164.000
	138.000
	174.000
	138.000
500.000	325.000
	442.000
	396.000
	369.000
1250.000	938.000
	952.000
	764.000
	902.000

TABLE II-132. NG on Transite - Analysis of Target-
Found Concentration Points

NITROGLYCERINE (NG)
TRANSITE SURFACE

ANALYSIS OF 16 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 531.25 SD= 450.693909433

FOUND CONC
MEAN= 377.421875 SD= 328.83558048

NO. RUNS 1 TOTAL X-Y ALL RUNS 16 NO. CONCENTR 16
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -3.90769230769
SLOPE= 0.721561538462

USE FOR ACCURACY

R= 0.988954389276

MEAN SQR DEV OF POINTS FROM REGRESSION= 2545.27895262
ST ERROR EST= 50.4507577011

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 14

TWO TAIL P LEVEL IS .1

t= 1.7613101065

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

Y(C)= 89.5956595952

X(D)= 260.722616729

NITROGLYCERINE (NG)
TRANSITE SURFACE
FOUND CONC

952.00




Figure 11-130. NG on Transite - Graph of Target-Found Concentration Points

TABLE II-133. NG on Transite - Inaccuracy and Imprecision Data

NITROGLYCERIN (NG)

TRANSITE SURFACE

STATISTICAL DATA USED TO DETERMINE PERCENT

INACCURACY AND IMPRECISION

Mn Target Con ug/10 sq cm	Mn Found Conc ug/10 sq cm	Standard Deviation	Mean Pct Inaccuracy	Imprecision
125.000	84.188	6.450	-32.650	7.661
250.000	153.500	18.358	-38.600	11.959
500.000	383.000	49.024	-23.400	12.800
1250.000	889.000	85.953	-28.880	9.669
Means		39.946	-30.883	18.522




Figure II-131. NG on Transite Graph of Inaccuracy

NITROGLYCERIN (NG)
TRANSITE SURFACE
IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-132. NG on Transite - Graph of Imprecision

TABLE II-134. PETN on Metal - Target vs. Found Concentrations

PENTERYTHRITE TETRAHITRATE (PETN)

PENTERYTHRITE TETRAHITRATE METAL SURFACE	
TARGET CONC. UG/ 10 SQ CM	US FOUND CONC Found Conc UG/ 10 SQ CM
50.000	30.200
50.000	32.400
50.000	47.400
50.000	47.800
100.000	78.000
100.000	73.500
100.000	90.800
100.000	79.600
200.000	199.000
200.000	159.000
200.000	187.000
200.000	189.000
500.000	450.000
500.000	406.000

Table 11-134. PETN on Metal - Target vs. Found Concentrations (Continued)

PENTAERYTHRITE TETRANITRATE (PETN)

METAL SURFACE	
TARGET CONC. UG/ 10 SQ CM	US FOUND CONC Found Conc UG/ 10 SQ CM
349.000	308.000
349.000	363.000
1000.000	902.000
1000.000	859.000
697.000	612.000
1388.000	1150.000

TABLE II-135. PETN on Metal - Analysis of Target-
Found Concentration Points

PENTAERYTHRITE TETRANITRATE (PETN)
METAL SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 359.15 SD= 382.937438205

FOUND CONC
MEAN= 313.185 SD= 328.387893672

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 5.97887894415

SLOPE= 0.855369959782

USE FOR ACCURACY

R= 0.997458150645

MEAN SQR DEV OF POINTS FROM REGRESSION= 577.940153813

ST ERROR EST= 24.0403858915

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

THO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 49.6274158122

x(d)= 101.531895278

PENTAERYTHRITE TETRANITRATE (PETN)
METAL SURFACE
FOUND CONC




Figure 11-133. PETN on Metal - Graph of Target-Found Concentration Points

TABLE II-136 PETN on Concrete - Target vs. Found Concentrations

(PETN)

PENTAERYTHRITE TETRANITRATE CONCRETE SURFACE	
TARGET CONC. UG/10 SQ CM	VS. FOUND CONC. Found Conc UG/10 SQ CM
50.000	30.600
50.000	44.400
50.000	45.800
50.000	44.200
100.000	97.000
100.000	86.100
100.000	84.700
100.000	86.900
200.000	167.000
200.000	194.000
200.000	181.000
200.000	182.000
349.000	323.000
349.000	292.000

Table II-136. PETN on Concrete - Target vs. Found Concentrations (Continued)

PENTAERYTHRITE TETRANITRATE (PETN)

CONCRETE SURFACE	
TARGET CONC. UG/ 10 SQ CM	FOUND CONC. UG/ 10 SQ CM
349.000	300.000
349.000	304.000
697.000	583.000
697.000	610.000
697.000	603.000
1388.000	1140.000

TABLE II-137 PEIN on Concrete - Analysis of
Target-Found Concentration Points

PENTAERYTHRITE TETRANITRATE (PETN)
CONCRETE SURFACE
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 313.75 SD= 333.434578922

FOUND CONC
MEAN= 270.335 SD= 276.82953766

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 10.0941892935

SLOPE= 0.829452783128

USE FOR ACCURACY

R= 0.999056104403

MEAN SQR DEV OF POINTS FROM REGRESSION= 152.635268141

ST ERROR EST= 12.3545646682

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F. = 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 32.5286725913

x(d)= 53.9135179075

PENTAERYTHRITE TETRANITRATE (PETN)
 CONCRETE SURFACE
 FOUND CONC

1149.00




Figure II-134. PETN on Concrete - Graph of Target-
 Found Concentration Points

TABLE II-138 PETN on Brick - Target vs. Found Concentrations

(PETN)

PENTAERYTHRITE TETRANITRATE	
BRICK SURFACE	
TARGET CONC. UG/ 10 SQ CM	VS FOUND CONC. Found Conc UG/ 10 SQ CM
50.000	39.700
50.000	25.000
50.000	26.400
50.000	22.200
100.000	70.400
100.000	47.200
100.000	52.800
100.000	47.200
200.000	112.000
200.000	100.000
200.000	97.100
200.000	101.000
500.000	364.000
349.000	184.000

Table II-138. PETN on Brick - Target vs. Found Concentrations (Continued)

PENTAERYTHRITE TETRANITRATE (PETN)

BRICK SURFACE	
TARGET CONC. UG/ 10 SQ CM	US FOUND CONC. Found Conc UG/ 10 SQ CM
349.000	168.000
349.000	180.000
1000.000	677.000
697.000	425.000
697.000	357.000
697.000	345.000

TABLE II-139 PETN on Brick - Analysis of Target-
Found Concentration Point

PENTAERYTHRITE TETRANITRATE (PETN)
BRICK SURFACE
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 301.9 SD= 277.200079745

FOUND CONC
MEAN= 172.05 SD= 174.438034542

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -14.2223284125
SLOPE= 0.617000094112
USE FOR ACCURACY
R= 0.980476968454
MEAN SQR DEV OF POINTS FROM REGRESSION= 1241.08232912
ST ERROR EST= 35.2403537031
USE FOR PRECISION
T FOR CONFIDENCE BAND
D.F.= 18
TWO TAIL P LEVEL IS .1
t= 1.73406096408
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C
Y(C)= 50.2303186945
X(D)= 206.251664127

PENTAERYTHRITE TETRANITRATE (PETN)
BRICK SURFACE
FOUND CONC

677.00




Figure 11-135. PETN on Brick - Graph of Target-Found Concentration Points

TABLE 11-140 PETN on Transite - Target vs. Found Concentrations

PENTERYTHRITE TETRANITRATE (PETN)

TRANSITE SURFACE TARGET CONC. US FOUND CONC. Target Conc Found Conc UG/ 10 SQ CM UG/ 10 SQ CM	
50.000	30.400
50.000	38.600
50.000	38.100
50.000	48.800
100.000	61.300
100.000	88.900
100.000	77.200
100.000	88.900
200.000	110.000
200.000	167.000
200.000	152.000
200.000	164.000
500.000	360.000
349.000	239.000

Table II-140. PETN on Transit - Target vs. Found Concentrations (Continued)

PENTAERYTHRITE TETRANITRATE (PETN)

TRANSIT SURFACE TARGET CONC. VS. FOUND CONC.	
Target Conc UG/ 10 SQ CM	Found Conc UG/ 10 SQ CM
349.000	288.000
349.000	299.000
1000.000	801.000
697.000	525.000
1000.000	804.000
1388.000	1158.000

TABLE 11-141 PETN on Transite - Analysis of Target-
Found Concentration Points

PENTAERYTHRITE TETRANITRATE (PETN)
TRANSITE SURFACE
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC
MEAN= 351.6 SD= 381.500141409

FOUND CONC
MEAN= 276.96 SD= 312.148352719

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -10.0036586066

SLOPE= 0.816165104114

USE FOR ACCURACY

R= 0.997497183375

MEAN SQR DEV OF POINTS FROM REGRESSION= 514.183809356

ST ERROR EST= 22.6756214767

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 31.1370879006

x(d)= 100.308711628

PENTAERYTHRITE TETRANITRATE (PETN)
 TRANSITE SURFACE
 FOUND CONC




Figure II-136. PETN on Transite - Graph of Target-Found Concentration Points